EN

SDB HbA1c Control
SDB HbA1c Control
PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

KIT CONTENTS

- SDB HbA1c Control-Level 1 : 10 tablets in the bottle (Transparent tube)
- SDB HbA1c Control-Level 2 : 10 tablets in the bottle (Red colored tube)
- Instructions for use

MATERIALS REQUIRED BUT NOT PROVIDED

- STANDARD F Analyzer
- STANDARD F HbA1c Test kit
- STANDARD A1Care Analyzer
- STANDARD A1Care Test Kit
- Multicare Analyzer
- Multicare HbA1c Test kit

CONTROL PREPARATION AND TEST PROCEDURE

This product should be treated the same as patient specimens and run in accordance with instruction accompanying the instrument and kit being used. It is recommended that low level and high level controls to be run:

- Once for each new lot
- Once for each untrained operator
- As required by test procedures in this instructions and in accordance with local, state and federal regulations of accreditation requirements

Preparation

- Expose STANDARD F HbA1c or STANDARD A1Care or Multicare HbA1c test kit and HbA1c control at room temperature (15-30°C/59-86°F) at least 30 minutes prior to the test.
- Carefully read instructions for use.
- Check the expiration date of the control on the bottle and of the test device on the pouch. Do not use expired controls or test devices.
- Open the foil pouch, and check the components.
- STANDARD F HbA1c
- Foil pouch
- spoit(Blue)
- Test device
- Result window
- Specimen well
- Desiccant

Control preparation

- Open the level 1 control bottle and take a tube out.
- Pour the control tablet into an extraction buffer tube in the test kit.
- Slightly shake the buffer tube to dissolve the tablet well.
- Insert the edge of the spoit into the extraction buffer tube to mix buffer with latex tablet in the spoit.
- Carefully press and release the rubber at the top of the spoit for 7-8 times to dissolve the latex tablet in the spoit. Try to minimize air bubbles while mixing.
- Collect all the mixture with the spoit from the extraction buffer tube.
- Too much bubbles in the spoit may induce erroneous results.

Analysis of control

[A1Care Test System]

- Prepare the A1Care test kit and A1Care analyzer. Place an A1Care analyzer on a level surface.
- Check the expiry date of the test kit. Do not use expired product.
- Insert the test panel into slot of the analyzer. When the test panel reaches the correct position, the test panel is automatically inserted deeper to the analyzer. If the test panel does not have any problem (damaged or defective), the blood and START button icons appear on the screen.
- Apply all control mixture into the test panel.
- After applying the control mixture, immediately press the 'TEST START' button.
- The analyzer will automatically display the test result after 3 minutes.
- Check the test result whether it is within the given range on the label of the control bottle. If the result is out of range, repeat the test. If the re-test value is also unacceptable, please contact the customer service center.
- Perform the same test with a Level 2 control tablet to see the performance of high value.

INTERPRETATION OF TEST RESULTS

HbA1c Level 1 Control		
Result	Interpretation	Follow up
Pass	Result within Level 1 control range	-
Fail	Result out of Level 1 control range	Retest
Invalid		Retest

HbA1c Level 2 Control		
Result	Interpretation	Follow up
Pass	Result within Level 2 control range	-
Fail	Result out of Level 2 control range	Retest
Invalid		Retest

[MultiCare Test System]

- Prepare the MultiCare HbA1c Test Kit and MultiCare Analyzer. Place a MultiCare Analyzer on a level surface.
- Insert a code chip into the analyzer and then insert a test panel with the 'HbA1c' inscription facing upwards. When the test panel reaches the correct position, the test panel is automatically inserted deeper into the analyzer. If the test panel does not have any problem (damaged or defective), the blood and START button icons appear on the screen.
- Check the given range on the each label of the control bottle.
- Use new test device and new control for retest.

EXPLANATION AND SUMMARY

Introduction

SDB HbA1c Control is designed to monitor the performance of STANDARD™ A1Care, MultiCare™ or STANDARD™ F System for HbA1c testing. To maintain the testing quality, at least two levels of appropriate control reagents should be tested on a regular basis. Quality control test should be performed in accordance with regulatory requirements and each laboratory's standard procedure.

Test principle

Two different levels of control help monitoring the performance of STANDARD™ A1Care, MultiCare™ or STANDARD™ F System for HbA1c testing. The range of each control is indicated on the bottle with a barcode which contains necessary information to process the quality control assessment. Both level 1 & level 2 test should be separately performed in order and results of tests should be within the given ranges.

ACTIVE INGREDIENTS OF MAIN COMPONENT

This product is composed of human hemoglobin protein with other preservatives and stabilizers. The range of control level is indicated on the label of each bottle.

- Level 1 control tablet - Human Hemoglobin 70 µg
- Level 2 control tablet - Glycated human hemoglobin 20 µg

KIT STORAGE AND STABILITY

Store the SDB HbA1c Control at 2-30°C/36-86°F. Do not freeze the Control material. The SDB HbA1c Control is stable until expiration date printed on the outer box.

WARNING AND PRECAUTION

- If there is evidence of microbial contamination in the reconstituted control, discard the control.
- Wear appropriate clothing and gloves when handling reagents.
- Clean any spillage by immediately and thoroughly wipe up with a suitable disinfectant such as 1% sodium hypochlorite solution.
- Handle and dispose of all controls and materials used in testing as though they contain infectious agents.
- Dispose of any discarded materials in accordance with the requirements of your local waste management authorities. In the event of damage to packaging, contact the local distributor of SD BIOSENSOR, Inc.
- The control bottle itself is a dehumidifier. Cap the control bottle tightly to maintain the quality.

LIMITATION OF TEST

- This product is provided for quality assurance purposes and must not be used for calibration or as primary reference preparations in any test procedure.
- Adverse storage conditions or use of outdated reagents may produce erroneous results.
- This product should not be used past the expiration date.
- Alterations in physical appearance may indicate instability or deterioration of this product. If there is evidence of microbial contamination in this product, discard it.

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