
**STANDARD Q COVID-19 Ag Test**

STANDARD™ Q COVID-19 Ag Test

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE PERFORMING THE TEST


SD BIOSENSOR

**KIT CONTENTS**

STANDARD COVID-19 Ag Positive Control swab  
STANDARD COVID-19 Ag Negative Control swab  
Instructions for use


**SPECIMEN COLLECTION AND PREPARATION****Specimen preparation (Nasal swab)**

1. Tilt patient's head back slightly.
  2. While rotating the swab, insert swab less than one inch (about 2cm) into nostril until resistance is met at turbinate.
  3. Rotate the swab 4 times against nasal wall.
  4. Repeat in other nostril using the same swab.
- Only one swab is used per test. Rub both nostrils with the same swab. Using proper pressure, rub the swab against all of the inside nasal walls of your nostril. Make at least 4 big circles along the nasal wall. Do not just spin the swab.


1. Insert the swab into an extraction buffer tube. While squeezing the buffer tube, stir the swab for 5 times.
2. Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.
3. Press the nozzle cap tightly onto the tube.
4. Specimen should be tested as soon as possible after collection.

1. Specimens may be stored at room temperature (15 – 25°C or 2-8°C/ 36-64°F for up to 4 hours prior to testing.
- Without the tube squeezing process, improper results may occur due to the large amount of buffer absorbed by the swab.
- If the specimen storage condition is out of instructions as below, do not use.
- The Nasal swab is stored in extraction buffer for more than 4 hours at 5±3°C or 20±5°C.
- Freezing and thawing of Nasal swab is more than 1 cycle.

**PREPARATION AND TEST PROCEDURE****Preparation**

1. Carefully read instructions for using STANDARD Q COVID-19 Ag Test.
2. Check the expiry date at the back of the foil pouch. Do not use the kit, if expiry date has passed.
3. Check the test device and the desiccant pack in the foil pouch.

**Test Procedure**

1. Apply 4 drops of extracted specimen to the specimen well of the test device.
2. Read the test result in 15-30 minutes.

- Place the test device on a flat surface.
- Dispense the specimen at 90 degree angle to allow for free falling drops and avoid bubbles.
- Do not read test results after 30 minutes. It may give false results.

**PERFORMANCE CHARACTERISTICS****Clinical evaluation**

Clinical performance of the STANDARD Q COVID-19 Ag Test was evaluated using nasal swab samples from 696 subjects in a prospective study at a clinical center in Germany. The study cohort included adults at high risk for SARS-CoV-2 infection according to clinical suspicion. 311 subjects underwent nasal sampling performed by healthcare professionals and 385 subjects followed instructions to obtain a nasal swab sample by themselves. Self-collection was performed under the supervision of healthcare workers without interference or assistance. Test procedures and result reading were always performed by healthcare professionals. PCR tests (Roche cobas® SARS-CoV-2 and TibMolbiol SARS-CoV-2 E-gene assay) using combined nasopharyngeal/oropharyngeal swab samples were used as the comparator methods. Nasal sampling always preceded the combined NP/OP sampling.

**Test sensitivity & specificity**

The following tables summarize the patient and performance characteristics of the STANDARD Q COVID-19 Ag Test. The relative sensitivity was 89.6% (CI: 30% to 95% CI: 79.7% - 95.7%) for professionally collected samples, and 89.1% (CI value < 30% to 95% CI: 78.8% - 95.5%) for self-collected samples. For patients for whom days post symptom onset was known, and was 0-5 days, the relative sensitivity in comparison to rPCR was 86.7% (95% CI: 75.4% - 94.1%) for professionally collected nasal samples and 88.9% (95% CI: 77.4% - 95.8%) for self-collected nasal samples. The relative specificity in comparison to rPCR was 99.1% (95% CI: 96.9% - 99.9%) for professionally collected nasal samples and 99.0% (95% CI: 97.2% - 99.8%) for self-collected nasal samples. In total, nasal swab samples from 150 PCR-positive and 546 PCR-negative individuals were evaluated using the STANDARD Q COVID-19 Ag Test. The relative sensitivity and relative specificity were 82.7% (95% CI: 75.6% - 88.4%) and 99.1% (95% CI: 97.3% - 99.7%), respectively.

Relative sensitivity, % (95% CI), N	Professional collection	Self-collection
C <sub>1</sub> ™ = 24	97.7% (CI: 88.0% - 99.9%), 44	97.9% (CI: 88.7% - 99.9%), 47
C <sub>2</sub> ™ = 27	93.1% (CI: 83.3% - 98.1%), 58	94.1% (CI: 85.4% - 98.9%), 57
C <sub>3</sub> ™ = 30	89.6% (CI: 79.7% - 95.7%), 67	89.1% (CI: 78.8% - 95.5%), 64
C <sub>4</sub> ™ = 33	87.1% (CI: 77.0% - 93.9%), 70	84.5% (CI: 74.0% - 92.0%), 71
All C values	83.1% (CI: 72.9% - 90.7%), 77	82.2% (CI: 71.5% - 90.2%), 73

depending upon the amount of SARS-CoV-2 antigen present in the specimen. If SARS-CoV-2 antigens are not present in the specimen, then no color appears in the test line. The control line is used for procedural control, and should always appear if the test procedure is performed properly and the test reagents of the control line are working.

**Kit contents (09COV33D)**


- 1 purple colored band will appear in the top section of the result window. This band is control line (C).
- A purple colored band will appear in the lower section of the result window. This band is test line of SARS-CoV-2 antigen (T).
- Even if the control line is faint, or the test line isn't uniform, the test should be considered to be performed properly and the test result should be interpreted as a positive result.

\* The presence of any line no matter how faint the result is considered positive. Positive results should be considered in conjunction with the clinical history and other data available.

**CONTROL PREPARATION AND TEST PROCEDURE****Positive/Negative control****Preparation**

1. Bring the kit contents and the specimens to room temperature before testing. Do not re-use the test kit.
2. Do not use the test kit if the pouch is damaged or the seal is broken.
3. Do not use the extraction buffer tube of another lot.
4. Do not smoke, drink or eat while handling specimens.
5. Wear personal protective equipment, such as gloves and lab coats when handling kit reagents. Wash hands thoroughly after the tests are done.
6. Clean up spills quickly using an appropriate disinfectant.
7. Handle all specimens as if they contain infectious agents.
8. Observe established precautions against microbiological hazards throughout testing procedures.
9. Dispose of all specimens and materials used to perform the test as biohazard waste. Laboratory chemical and biohazard wastes must be biohazard and discarded in accordance with all local, state, and national regulations. Desiccant in foil pouch is to absorb moisture and keep humidity from affecting products. If the moisture indicating desiccant beads change from yellow to green, the test device in the pouch should be discarded.
10. This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

- Warning: H317 May cause an allergic skin reaction. H412 Harmful to aquatic life with long lasting effects. H319 Causes serious eye irritation.
- Prevention: P261 Avoid breathing dust/fume/gas/mist/vapours/spray. P273 Avoid release to the environment. P280 Wear eye protection/face protection.

**Test Procedure**

1. Put the positive or negative control swab into an extraction buffer tube. While squeezing the buffer tube, stir the swab more than 5 times.
2. Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.
3. Press the nozzle cap tightly onto the tube.

1. Apply 4 drops of extraction buffer to the specimen well.
2. Read the test result in 15-30 minutes.

- Place the test device on a flat surface.
- Dispense the specimen at 90 degree angle to allow for free falling drops and avoid bubbles.
- Do not read test results after 30 minutes. It may give false results.

**INTERPRETATION OF CONTROL TEST RESULT**

STANDARD COVID-19 Ag Control Positive swab: Positive	Result	Interpretation	Follow up
Test (T) Line Positive	PASS	-	-
Test (T) Line Negative	FAIL	Retest	-
No Control (C) Line	Invalid	-	Retest

STANDARD COVID-19 Ag Negative Control swab: Negative	Result	Interpretation	Follow up
Test (T) Line Negative	PASS	-	-
Test (T) Line Positive	FAIL	Retest	-
No Control (C) Line	Invalid	-	Retest

**EXPLANATION AND SUMMARY****Introduction**

Coronavirus is a single-stranded positive-sense RNA virus with an envelope of about 80 to 120 nm in diameter. Its genetic material is the largest of all RNA viruses and is an important pathogen of many domestic animals, pets, and human diseases. It can cause a variety of acute and chronic diseases. Common signs of a person infected with a coronavirus include respiratory symptoms, fever, cough, shortness of breath, and dyspnea. In more severe cases, infection can cause pneumonia, severe acute respiratory syndrome, kidney failure, and even death. The 2019 new coronavirus, or SARS-CoV-2 (COVID-19), was discovered because of Wuhan Viral Pneumonia cases in 2019, and was named by the World Health Organization on January 12, 2020, confirming that it can cause colds and the Middle East Respiratory Syndrome (MERS) and more serious diseases such as Severe Acute Respiratory Syndrome (SARS). These kits are helpful for the auxiliary diagnosis of coronavirus infection. The test results are for clinical reference only and cannot be used as a basis for confirming or excluding cases alone.

**Intended use**

STANDARD Q COVID-19 Ag Test is a rapid chromatographic immunoassay for the qualitative detection of specific nucleocapsid protein antigen from SARS-CoV-2 present in human nasal or nasopharyngeal specimens. This product is intended for healthcare professionals at the clinical setup and point of care sites, as an aid to early diagnosis of SARS-CoV-2 infection in patient with clinical symptoms of SARS-CoV-2 infection. It provides only an initial screening test result. This product is strictly for medical professional use only and not intended for personal use. The administration of the test and the interpretation of the results should be done by a trained health professional. The result of this test should not be the sole basis for the diagnosis; confirmatory testing is required.

**Test principle**

STANDARD Q COVID-19 Ag test device has two pre-coated lines, "C" Control line, "T" Test line on the surface of the nitrocellulose membrane. Both the control line and test line in the result window are not visible before applying any specimens. Mouse monoclonal anti-SARS-CoV-2 antibody is coated on the test line region and mouse monoclonal anti-Chicken IgY antibody is coated on the control line region. Mouse monoclonal anti-SARS-CoV-2 antibody conjugated with color particles are used as detectors for SARS-CoV-2 antigenic. During the test, SARS-CoV-2 antigen in the specimen interact with the monoclonal anti-SARS-CoV-2 antibody conjugated with color particles making antigen-antibody color particle complex. This complex migrates on the membrane via capillary action until the test line, where it will be captured by the mouse monoclonal anti-SARS-CoV-2 antibody. A colored test line would be visible in the result window if SARS-CoV-2 antigens are present in the specimen. The intensity of colored test line will vary

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**MATERIALS REQUIRED BUT NOT PROVIDED**

1. Personal protective Equipment per local recommendations (i.e. gown/lab coat, face mask, face shield/eye goggles and gloves)
2. Timer
3. Biohazard container

**KIT STORAGE AND STABILITY**

Store the kit at 2-30°C / 36-86°F out of direct sunlight. Kit materials are stable until the expiration date printed on the outer box. Do not freeze the kit.

**WARNINGS AND PRECAUTIONS**

1. Bring the kit contents and the specimens to room temperature before testing. Do not re-use the test kit.
2. Do not use the test kit if the pouch is damaged or the seal is broken.
3. Do not use the extraction buffer tube of another lot.
4. Do not smoke, drink or eat while handling specimens.
5. Wear personal protective equipment, such as gloves and lab coats when handling kit reagents. Wash hands thoroughly after the tests are done.
6. Clean up spills quickly using an appropriate disinfectant.
7. Handle all specimens as if they contain infectious agents.
8. Observe established precautions against microbiological hazards throughout testing procedures.
9. Dispose of all specimens and materials used to perform the test as biohazard waste. Laboratory chemical and biohazard wastes must be biohazard and discarded in accordance with all local, state, and national regulations. Desiccant in foil pouch is to absorb moisture and keep humidity from affecting products. If the moisture indicating desiccant beads change from yellow to green, the test device in the pouch should be discarded.
10. This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

- Warning: H317 May cause an allergic skin reaction. H412 Harmful to aquatic life with long lasting effects. H319 Causes serious eye irritation.
- Prevention: P261 Avoid breathing dust/fume/gas/mist/vapours/spray. P273 Avoid release to the environment. P280 Wear eye protection/face protection.

- For customers in the European Economic Area: Contains SVHC: octyl/nonylphenol ethoxylates.
- For use as part of an IVD method and under controlled conditions only – acc. to Art. 56.3 and 3.23 REACH Regulation.

**LIMITATION OF TESTS**


1. The test procedure, precautions and interpretation of results for this test must be followed strictly when testing.
2. The test should be used for the detection of SARS-CoV-2 antigen in human nasal swab specimens only, other specimen types have not been validated.
3. This test can not be used for quantifying SARS-CoV-2 antigen concentration. Failure to follow the test procedure and interpretation of test results may adversely affect test performance and/or produce invalid results.
4. The test result must always be evaluated with other data available to the physician.
5. A negative result may occur if the concentration of antigen in a specimen is below the detection limit of the test or if the specimen was collected or transported improperly, therefore a negative test result does not eliminate the possibility of SARS-CoV-2 infection, and should be confirmed by viral culture or molecular assay.
6. Positive test results do not rule out co-infections with other pathogens.

ES **REF Q-NCOV-016****STANDARD Q COVID-19 Ag Test**

STANDARD™ Q COVID-19 Ag Test

PLEASE READ THE INSTRUCTIONS CAREFULLY BEFORE REALIZING THE PUEBA.

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Summary of sample characteristics and performances				
	Overall	HCP-collection	Self-collection	
N	696	311	385	
Asymptomatic, n/N (%)	20/696 (2.9%)	7/311 (2.3%)	13/385 (3.4%)	
Symptomatic, n/N (%)	676/696 (97.1%)	304/311 (97.7%)	372/385 (96.6%)	
DP50, median (range)	3 (0 - 27)	3 (0 - 15)	4 (0 - 27)	
PCR positive, n/N (%)	150/696 (21.6%)	77/311 (24.8%)	73/385 (19.0%)	
PCR positive symptomatic, n/N (%)	147/580 (98.0%)	75/77 (97.4%)	72/73 (98.6%)	
PCR positive asymptomatic, n/N (%)	3/150 (2.0%)	2/77 (2.6%)	1/73 (1.4%)	
PCR negative, n/N (%)	546/696 (78.4%)	234/311 (75.2%)	312/385 (81.0%)	
Relative sensitivity, % (95% CI), N		Combined DP/NP		
C <sub>1</sub> ™ = 24	97.7% (CI: 88.0% - 99.9%), 44		97.9% (CI: 88.7% - 99.9%), 47	
C <sub>2</sub> ™ = 27	93.1% (CI: 83.3% - 98.1%), 58		94.1% (CI: 85.4% - 98.9%), 57	
C <sub>3</sub> ™ = 30	89.6% (CI: 79.7% - 95.7%), 67		89.1% (CI: 78.8% - 95.5%), 64	
C <sub>4</sub> ™ = 33	87.1% (CI: 77.0% - 93.9%), 70		84.5% (CI: 74.0% - 92.0%), 71	
All C values	83.1% (CI: 72.9% - 90.7%), 77		82.2% (CI: 71.5% - 90.2%), 73	
a) for samples run on cobas the Target E (Gene) Ct values were used.				
Relative specificity, % (95% CI), N		Professional collection	Self-collection	
All Ct values	99.1% (CI: 96.9% - 99.9%), 234		99.0% (CI: 97.2% - 99.8%), 312	

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© SD BIOSENSORS