STANDARD F HbA1c



» Background

STANDARD F HbA1c is for quantitative measurement of glycated hemoglobin (HbA1c) in human capillary or venous whole blood. This test is to monitor glycemic control in people with diabetes.

» Product Specification









	The state of the s		
STANDARD F HbA1c			
Intended use	Quantitative measurement of glycated hemoglobin (HbA1c)		
Specimen type	Whole blood		
Specimen volume	5 μl		
Measurement range	4 – 15 % [NGSP], 20 – 140 mm/mol [IFCC]		
Reference range	≤ 5.6% (Normal) 5.7 – 6.4% (Prediabetes) ≥ 6.5% (Diabetes) 7% (ADA target for diabetes patients)		
Testing time	3 mins		
No Interference	Common Hb variants (HbS, HbE, HbC, HbD, HbF)		
Storage condition	2 – 30°C / 36 - 86°F		

» Test Procedure



Collect 5 μ I of blood sample with the SpoitTM (Blue).




Mix the sample with the extraction buffer and collect all the mixture with the Spoit TM (Blue).



Apply the sample mixture into the sample well of the test device.

» Method Comparison



Reference method vs STANDARD F HbA1c		
Correlation with HPLC	y=0.9628x - 0.2768; R ² =0.9906, n=80	
CV%	< 3%	
Differ(%)	within 6% (NGSP criteria)	

» Ordering Information

Category	Product	Pack Size	CAT No.
Chronic disease —	STANDARD F HbA1c	20 Tests	10A1C10B
	SDB HbA1c Control	Lv1 x 10 / Lv2 x 10	03ACS10