
STANDARD F β-hCG FIA

STANDARD™ β-hCG FIA
PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

MATERIALS REQUIRED BUT NOT PROVIDED

- STANDARD F Analyzer

SPECIMEN COLLECTION AND PREPARATION

Whole blood

- Collect the venous whole blood into the commercially available EDTA tube by venipuncture.
- It is recommended that collected venous whole blood specimens are used immediately. If venous whole blood in an anticoagulant tube is stored at room temperature and refrigerator condition (2-8°C/36-46°F), the specimen can be used for testing within 8 hours after collection.
- Do not use hemolyzed or frozen blood specimens.


Serum

- Collect the whole blood into the commercially available plain tube. NOT containing anti-coagulants such as heparin, EDTA or sodium citrate by venipuncture and leave to settle for 30 minutes for blood coagulation and then centrifuge blood to get serum specimen of supernatant.
- Serum specimen may be stored at room temperature for up to 8 hours after collection and in refrigerator conditions (2-8°C/36-46°F) for up to 3 days prior to testing.
- It should be brought to room temperature prior to use.

TEST PROCEDURE

Preparation

- Allow kit components and collected specimen to room temperature at least 30 minutes before starting the test.
- Carefully read instructions for the STANDARD F β-hCG FIA.
- Check the expiry date at the back of the foil pouch. Use another lot if expiry date has passed.
- Open the foil pouch and check the test device in the foil pouch.

Analysis of specimen


STANDARD TEST mode STANDARD F100, F200 and F2400 Analyzers

- Take the test device out of the foil pouch and place it on a flat and dry surface. Write patient information on the label of test device.

STANDARD F2400 analyzer	Workplace	Run Test	Scan or type patient ID and/or operator ID
STANDARD F100 and F200 analyzer	Standard Test mode	Insert patient ID and / or operator ID on the analyzer	

- Insert the test device to the test slot of the analyzer. When inserting the test device to the analyzer, the analyzer read the barcode data, and check the test device is valid.

- Collect 50µl of specimen with a fixed volume pipette (50µl) - [①] labeled.

PERFORMANCE CHARACTERISTICS

Analytical Performance

1. Accuracy (Method Comparison)

Results comparing the STANDARD™ β-hCG FIA with an automated chemiluminescent method are shown below.

1) Serum

Regression Analysis	
Slope	1.0161
Y-Intercept	-6.6452
R	0.9973
R ²	0.9947

2) Whole blood

Regression Analysis	
Slope	0.9889
Y-Intercept	2.2289
R	0.9979
R ²	0.9958

3. Precision

For precision evaluation, the within-run using 3 levels and the Day-to-Day using 3 levels for 20 days were done at 3 different testing sites. The acceptance criterion were within 10%(CV) for serum and 15%(CV) for whole blood.

- Do not use any artificial materials.
- Place the analyzer on a flat surface when in use.
- Wash your hands in warm and soapy water. Rinse well and dry completely before testing.
- Discard the used test kit according to the proper method.
- Desiccant in foil pouch is to absorb moisture and keep humidity from affecting products. If the moisture indicating desiccant beads change from yellow to green, the test device in the pouch should be discarded.
- Use a fixed volume pipette (50µl) - [①] labeled on the pouch for specimen extraction purpose only. Do not use it as a specimen dispenser.
- Use a disposable dropper (100µl) - [②] labeled on the pouch for mixing the specimen, and for dispensing the specimen mixture into the test device.
- Check the expiration date printed at the pouch or package.
- Use the volume (150µl) of extraction buffer.
- Use the STANDARD F β-hCG FIA at room temperature.
- All kit components are must be at room temperature 30 minutes before running the assay.
- Do not write on the barcode or damage the barcode of the test device.

LIMITATION OF TEST

- The test should be used for the detection of β-hCG levels in human whole blood and serum specimens.
- Failure to follow the test procedure and interpretation of test results may adversely affect test performance and/or produce invalid results.
- Invalid results may occur if a poor-quality specimen is obtained.
- The test result must always be evaluated with other clinical data available to the physician.

QUALITY CONTROL

Calibration

The calibration set test of STANDARD F analyzer should be conducted according to the analyzer's manual.


When to use calibration set

- Before using the analyzer for the first time.
- When you drop the analyzer.
- Whenever you do not agree with the final result.
- When you want to check the performance of the analyzer and test device.

How to use calibration set

Calibration set test is a required function that ensures optimal performance by checking the internal analyzer optics and functions.

- Select the "Calibration" in main menu.
- The specific calibration set is included with the analyzer.
- Insert the CAL-1 for white calibration, CAL-2 for UV LED calibration, and CAL-3 for RGB LED calibration in sequence.

INTERPRETATION OF TEST RESULTS

Measuring Range 5 - 1,500 mIU/mL

Out of measuring range If below 5 mIU/mL, 'S 1' message will be displayed. If above 1,500 mIU/mL, '1,500 1' message will be displayed.

Early pregnancy diagnosis

Weeks since LMP	β-hCG Levels (mIU/mL)
3	5-50
4	5-426
5	18-7,340
6	1,080-56,500

The β-hCG reference ranges are provided for orientational purpose only.

Clinicians should use the test results in conjunction with the patient's other diagnostic findings and clinical signs, and interpret the concrete values in the context of the patient's clinical situation.

EXPLANATION AND SUMMARY

Introduction

Human chorionic gonadotropin (hCG) is a glycoprotein hormone, with structural similarity to the pituitary hormones follicle stimulating hormone (FSH), thyroid stimulating hormone (TSH), and luteinizing hormone (LH). All are glycoproteins consisting of two non-covalently bound dissimilar subunits, designated alpha and beta, with attached carbohydrate side chains. The alpha subunits of these glycoproteins are very similar. But the beta subunits differ, and confer immunological and biological specificity. Shortly after implantation of a fertilized ovum into the uterine wall, the trophoblast begins to produce hCG, which maintains steroid secretions of the corpus luteum until the placenta can do so. hCG can be detected after implantation; concentrations double approximately every 1.5-2 days for the first 6 weeks and then continue to rise until the end of the first trimester, gradually falling to a lower level for the remainder of the pregnancy. After delivery, hCG is usually undetectable several days postpartum. The hormone is an excellent marker for pregnancy. Healthy, non-pregnant individuals have low to undetectable hCG levels; however, hCG, originating from the pituitary gland, can be found at detectable levels in peri- and post-menopausal women. During pregnancy, unusually low or rapidly declining levels may indicate an abnormal condition such as an ectopic pregnancy or impending spontaneous abortion. STANDARD F β-hCG FIA, employing immunofluorescent detection system with STANDARD F analyzer, provides significantly fast, easy and accurate system to screen for pregnancy.

Intended use

STANDARD F β-hCG FIA is a fluorescence immunoassay for the quantitative determination of total β-hCG levels in human serum and whole blood using STANDARD F Analyzers manufactured by SD BIOSENSOR. This test is an *in vitro* diagnostic use and intended for use as an aid in the early detection of pregnancy.

Test Principle

When a human specimen is applied into the specimen well of the test device, the specimen migrates through the membrane from the specimen well. If β-hCG is present, it will be bound to by detector coupled to europium microparticle that migrates through the membrane. The europium microparticle complex will be captured by capture antibody on the test line where it is detected by STANDARD F analyzer. If β-hCG is not present, the europium microparticle will not be trapped by the capture antibody nor detected by STANDARD F Analyzer. The intensity of the fluorescence light generated on the membrane is scanned by the STANDARD F Analyzer. STANDARD F Analyzer can analyze the presence of the β-hCG in the clinical specimen by processing the results using pre-programmed algorithms and display the test result on the screen.

KIT STORAGE AND STABILITY

For 18 months when stored at 2-30°C/36-86°F and 10-90% RH.

WARNINGS AND PRECAUTIONS

- STANDARD F β-hCG FIA is for *in vitro* diagnostics use only.
- Carefully follow the instructions and procedures described in this Instructions for use before testing.
- STANDARD F β-hCG FIA should be used with STANDARD F Analyzer.
- STANDARD F β-hCG FIA should remain in its original sealed pouch until ready to use. Do not use the test device if pouch is damaged or the seal is broken.
- STANDARD F β-hCG FIA is single use only. Do not re-use it.
- Do not use hemolyzed specimens or frozen specimens.

- Do not use any artificial materials.
- Place the analyzer on a flat surface when in use.
- Wash your hands in warm and soapy water. Rinse well and dry completely before testing.
- Discard the used test kit according to the proper method.
- Desiccant in foil pouch is to absorb moisture and keep humidity from affecting products. If the moisture indicating desiccant beads change from yellow to green, the test device in the pouch should be discarded.
- Use a fixed volume pipette (50µl) - [①] labeled on the pouch for specimen extraction purpose only. Do not use it as a specimen dispenser.
- Use a disposable dropper (100µl) - [②] labeled on the pouch for mixing the specimen, and for dispensing the specimen mixture into the test device.
- Check the expiration date printed at the pouch or package.
- Use the volume (150µl) of extraction buffer.
- Use the STANDARD F β-hCG FIA at room temperature.
- All kit components are must be at room temperature 30 minutes before running the assay.
- Do not write on the barcode or damage the barcode of the test device.

LIMITATION OF TEST

- The test should be used for the detection of β-hCG levels in human whole blood and serum specimens.
- Failure to follow the test procedure and interpretation of test results may adversely affect test performance and/or produce invalid results.
- Invalid results may occur if a poor-quality specimen is obtained.
- The test result must always be evaluated with other clinical data available to the physician.

QUALITY CONTROL

Calibration

The calibration set test of STANDARD F analyzer should be conducted according to the analyzer's manual.


When to use calibration set

- Before using the analyzer for the first time.
- When you drop the analyzer.
- Whenever you do not agree with the final result.
- When you want to check the performance of the analyzer and test device.

How to use calibration set

Calibration set test is a required function that ensures optimal performance by checking the internal analyzer optics and functions.

- Select the "Calibration" in main menu.
- The specific calibration set is included with the analyzer.
- Insert the CAL-1 for white calibration, CAL-2 for UV LED calibration, and CAL-3 for RGB LED calibration in sequence.

INTERPRETATION OF TEST RESULTS

Measuring Range 5 - 1,500 mIU/mL

Out of measuring range If below 5 mIU/mL, 'S 1' message will be displayed. If above 1,500 mIU/mL, '1,500 1' message will be displayed.

Early pregnancy diagnosis

Weeks since LMP	β-hCG Levels (mIU/mL)
3	5-50
4	5-426
5	18-7,340
6	1,080-56,500

The β-hCG reference ranges are provided for orientational purpose only.

Clinicians should use the test results in conjunction with the patient's other diagnostic findings and clinical signs, and interpret the concrete values in the context of the patient's clinical situation.

EXPLANATION AND SUMMARY

Introduction

Human chorionic gonadotropin (hCG) is a glycoprotein hormone, with structural similarity to the pituitary hormones follicle stimulating hormone (FSH), thyroid stimulating hormone (TSH), and luteinizing hormone (LH). All are glycoproteins consisting of two non-covalently bound dissimilar subunits, designated alpha and beta, with attached carbohydrate side chains. The alpha subunits of these glycoproteins are very similar. But the beta subunits differ, and confer immunological and biological specificity. Shortly after implantation of a fertilized ovum into the uterine wall, the trophoblast begins to produce hCG, which maintains steroid secretions of the corpus luteum until the placenta can do so. hCG can be detected after implantation; concentrations double approximately every 1.5-2 days for the first 6 weeks and then continue to rise until the end of the first trimester, gradually falling to a lower level for the remainder of the pregnancy. After delivery, hCG is usually undetectable several days postpartum. The hormone is an excellent marker for pregnancy. Healthy, non-pregnant individuals have low to undetectable hCG levels; however, hCG, originating from the pituitary gland, can be found at detectable levels in peri- and post-menopausal women. During pregnancy, unusually low or rapidly declining levels may indicate an abnormal condition such as an ectopic pregnancy or impending spontaneous abortion. STANDARD F β-hCG FIA, employing immunofluorescent detection system with STANDARD F analyzer, provides significantly fast, easy and accurate system to screen for pregnancy.


Intended use

STANDARD F β-hCG FIA is a fluorescence immunoassay for the quantitative determination of total β-hCG levels in human serum and whole blood using STANDARD F Analyzers manufactured by SD BIOSENSOR. This test is an *in vitro* diagnostic use and intended for use as an aid in the early detection of pregnancy.

Test Principle

When a human specimen is applied into the specimen well of the test device, the specimen migrates through the membrane from the specimen well. If β-hCG is present, it will be bound to by detector coupled to europium microparticle that migrates through the membrane. The europium microparticle complex will be captured by capture antibody on the test line where it is detected by STANDARD F analyzer. If β-hCG is not present, the europium microparticle will not be trapped by the capture antibody nor detected by STANDARD F Analyzer. The intensity of the fluorescence light generated on the membrane is scanned by the STANDARD F Analyzer. STANDARD F Analyzer can analyze the presence of the β-hCG in the clinical specimen by processing the results using pre-programmed algorithms and display the test result on the screen.

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BIBLIOGRAPHY

- Sturgeon CM and McAllister EJ. Analysis of hCG: clinical applications and assay requirements. Ann Clin Biochem 1996; 35: 460-491.
- National Committee for Clinical Laboratory Standards. Chorionic gonadotropin Testing: Nomenclature, Reference Ranges, Assay Performance, and Clinical Application; Approved Guideline NCCLS Document 1/A10-A, December 1996.
- Kadar N, Caldwell BV, Romero R. A Method of Screening for Ectopic Pregnancy and Its Indications. Obstet Gynecol 1981; 58:162-6.
- Lenton EA, Neal LM, Sulaiman R. Plasma Concentrations of Human Chorionic Gonadotropin from the Time of Implantation until the Second Week of Pregnancy. Fertil Steril 1982; 37:773-8.
- Gronowski AM, et al. Use of Serum FSH to Identify Perimenopausal Women with Pituitary hCG. Clin Chem 2008; 54 (4): 652-656.

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