
STANDARD F COVID-19 IgM/IgG Combo FIA

STANDARD™ F COVID-19 IgM/IgG Combo FIA
PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

KIT CONTENTS

SPECIMEN COLLECTION AND PREPARATION

Serum
1. Collect the whole blood into the commercially available plain tube. NOT containing anti-coagulants such as heparin, EDTA or citrate, by venipuncture and allow to settle for 30 minutes for blood coagulation and then centrifuge blood to get serum specimen of supernatant.
2. If serum in the plain tube is stored in a refrigerator at 2-8°C/36-46°F, the specimen can be used for testing within 1 week after collection. Using the specimen in the long-term keeping more than 1 week can cause non-specific reaction. For prolonged storage, it should be stored at below -40°C/-40°F. They should be brought to room temperature prior to use.

Plasma
1. Collect the venous blood into the commercially available EDTA anti-coagulant tube by venipuncture and centrifuge blood to get plasma specimen.
2. If plasma in an anti-coagulant tube is stored in a refrigerator at 2-8°C/36-46°F, the specimen can be used for testing within 1 week after collection. Using the specimen in the long-term keeping more than 1 week can cause non-specific reaction. For prolonged storage, it should be stored at below -40°C/-40°F. They should be brought to room temperature prior to use.

Whole blood
[Capillary whole blood]
1. Capillary whole blood should be collected aseptically by fingertip.
2. Clean the area to be lanced with an alcohol swab.
3. Squeeze the end of the fingertip and pierce with a sterile lancet.
4. Using a capillary tube, collect the 20μl of capillary whole blood to the black line of the capillary tube.
5. The capillary whole blood must be tested immediately after collection.

[Venous whole blood]
1. Collect the venous whole blood into the commercially available anti-coagulant tube such as heparin.
2. If venous whole blood in an anti-coagulant tube is stored in a refrigerator at 2-8°C/36-46°F, the specimen can be used for testing within 1 day after collection.

3. Do not use hemolyzed blood specimens.


Use separate disposable materials for each specimen in order to avoid cross-contamination which can cause erroneous results.

PREPARATION AND TEST PROCEDURE

Preparation
1. Allow test device and collected specimen to room temperature (15-30°C / 59-86°F) prior to testing.
2. Carefully read instructions for using the STANDARD F COVID-19 IgM/IgG Combo FIA.
3. Take out the test device in foil pouch.

* Check the valid expiry date at the back of the foil pouch. Do not use if the expiry date has passed.

* The valid color of indicator in desiccant is yellow.


The violet-line on the membrane of unused test device will disappear after use.

Do not write on the barcode or damage the barcode of the test device.

■ Analysis of specimen

STANDARD TEST mode	STANDARD F100, F200 and F2400 analyzer
1. Prepare a STANDARD F Analyzer and select the 'Standard Test' mode according to the analyzer's manual.	STANDARD F100, F200 and F2400 analyzer 'Workplace' → 'Run Test' → Insert patient ID and / or operator ID on the analyzer 'Standard Test' mode → Insert patient ID and / or operator ID on the analyzer

2. Take the test device out of the foil pouch and place it on a flat and dry surface. Write patient information on the label of test device.

3. Insert the test device to the test slot of the analyzer. When inserting the test device to the analyzer, the analyzer will read the barcode data, and check the test device is valid.

4. Add the collected serum and plasma (10μl) or whole blood (20μl) to the specimen well of the test device.

5. Apply 3 drops of buffer to the specimen well in the test device.

6. After applying the buffer, immediately press the 'TEST START' button.

7. The analyzer will automatically display the test result in 15 minutes.

'READ ONLY' mode

1. Take the test device out of the foil pouch and place it on a flat and dry surface. Write information on the label of test device.

2. Add the collected serum and plasma (10μl) or whole blood (20μl) to the specimen well of the test device.

3. Apply 3 drops of buffer to the specimen well in the test device.

4. Incubate the test device for 15 minutes outside of the analyzer.

5. Prepare a STANDARD F Analyzer and select the 'Read Only' mode according to the analyzer's manual.

6. Insert the test device to the test slot of the analyzer. The analyzer automatically reads the barcode information of the test device and outputs the resulting values based on the test information.

7. Positive results should be considered in conjunction with the clinical history and other data available to the physician.

INTERPRETATION OF TEST RESULT

Result	COI (Cutoff index) value	Interpretation
Positive	COI ≥ 1.0	Positive result for Anti-SARS-CoV-2 IgM and/or IgG
Negative	COI < 1.0	Negative result for Anti-SARS-CoV-2 IgM and/or IgG
Invalid	Not show the COI value	Retest should be performed

Positive results should be considered in conjunction with the clinical history and other data available to the physician.

Do not write on the barcode or damage the barcode of the test device.

PERFORMANCE CHARACTERISTICS

Clinical evaluation
Performance characteristic for the STANDARD F COVID-19 IgM/IgG Combo FIA for rapid detection of anti-SARS-CoV-2 antibodies was established in retrospective, single institute, randomized, single-blinded study conducted at a trial site in KOREA during the 2020 SARS-CoV-2 pandemic situation. A total of 280 retrospective specimens were tested using the STANDARD F COVID-19 IgM/IgG Combo FIA with STANDARD F Analyzer. These specimens consisted of serum from PCR positive or negative confirmed patients. The performance of the STANDARD F COVID-19 IgM/IgG Combo FIA were compared to a commercialized molecular assay. Although the STANDARD F COVID-19 IgM/IgG Combo FIA allows to test for IgM and IgG separately, due to the differing inter-patient time response to the virus, any individual with positive result for the IgM or the IgG test should be read as a positive for anti-SARS-CoV-2 antibodies. The combined test result (positive for IgM and/or IgG or negative for IgM and/or IgG) was used to calculate the total test sensitivity and specificity.

[Test sensitivity and specificity]

The seroconversion time of IgM and IgG antibodies varies from person to person, but it was estimated to be around 7 days after onset of symptom. The STANDARD F COVID-19 IgM/IgG Combo FIA showed 100% (119/119) of sensitivity using specimens from patients 7 days after symptom onset (combined IgM+IgG). The STANDARD F COVID-19 IgM/IgG Combo FIA showed 95.33% (143/150) of specificity (combined IgM+IgG).

Table 2. Summary of the sensitivity and specificity of the STANDARD F COVID-19 IgM/IgG Combo FIA compared to PCR using specimens (COVID-19 IgM)		
PCR		
	Positive	Negative
STANDARD F COVID-19 IgM/IgG Combo FIA	120	3
Positive	120	3
Negative	10	147
Total	130	150
Sensitivity	92.31% (120/130, 95% CI: 86.31% - 96.25%)	
Specificity	98.00% (147/150, 95% CI: 94.27% - 99.59%)	

Table 4. Summary of the sensitivity of the STANDARD F COVID-19 IgM/IgG Combo FIA compared to PCR using specimens within less than 7 days after symptom onset (combined COVID-19 IgM+IgG)		
PCR		
	Positive	Negative
STANDARD F COVID-19 IgM/IgG Combo FIA	9	0
Positive	9	0
Negative	2	0
Total	11	0
Sensitivity	81.82% (9/11, 95% CI: 48.22% - 97.72%)	

Table 5. Summary of the sensitivity of the STANDARD F COVID-19 IgM/IgG Combo FIA compared to PCR using specimens 7 days after symptom onset (combined COVID-19 IgM+IgG)		
PCR		
	Positive	Negative
STANDARD F COVID-19 IgM/IgG Combo FIA	119	0
Positive	119	0
Negative	0	0
Total	119	0
Sensitivity	100% (119/119, 95% CI: 96.95% - 100%)	

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8. Positive test results do not rule out co-infections with other pathogens.
9. Negative test results are not intended to rule in other coronavirus infection except the SARS-CoV.
10. Children tend to shed virus for longer periods of time than adults, which may result in differences in sensitivity between adults and children.

QUALITY CONTROL

■ STANDARD F Analyzers calibration check

The calibration set test of STANDARD F Analyzers should be conducted according to the analyzer's manual.

[When to use calibration set]

1. Before using the analyzer and select the 'Standard Test' mode according to the analyzer's manual.
2. STANDARD F2400 analyzer
'Workplace' → 'Run Test' → Insert patient ID and / or operator ID on the analyzer
3. STANDARD F100 and F200 analyzer
'Standard Test' mode → Insert patient ID and / or operator ID on the analyzer
4. Take the test device out of the foil pouch and place it on a flat and dry surface. Write patient information on the label of test device.
5. Insert the test device to the test slot of the analyzer. When inserting the test device to the analyzer, the analyzer will read the barcode data, and check the test device is valid.
6. Add the collected serum and plasma (10μl) or whole blood (20μl) to the specimen well of the test device.
7. Add 3 drops of buffer to the specimen well in the test device.
8. After applying the buffer, immediately press the 'TEST START' button.
9. The analyzer will automatically display the test result in 15 minutes.

The STANDARD F Analyzer automatically calibrates and identifies the optical performance through measuring the membrane of the test device whenever the test is conducted in the 'Standard Test' mode. If the 'EEF' message displays on the screen, it means that the analyzer has a problem, so check with CAL devices. Contact the SD BIOSENSOR local distributor for further assistance.

[How to use calibration set]

Calibration set test is a mandatory function that ensures optimal performance by checking the internal analyzer optics and functions.

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[What is the calibration set purpose?]

The calibration set purpose is to help the auxiliary diagnosis of SARS-CoV-2 infection. The test results are for clinical reference only and cannot be used as a basis for confirming or excluding cases alone.

[Intended use]

STANDARD F COVID-19 IgM/IgG Combo FIA is the fluorescent immunoassay for the qualitative detection of specific antibodies to SARS-CoV-2 present in human serum, plasma and whole blood. STANDARD F COVID-19 IgM/IgG Combo FIA should be used with the STANDARD F Analyzers manufactured by SD BIOSENSOR.

This test is for in vitro diagnostic use and intended as an aid to diagnosis of SARS-CoV-2 infection in convalescent phase of patient with clinical symptoms with SARS-CoV-2 infection. It provides only an initial screening test result. More specific alternative diagnosis methods should be performed in order to obtain the confirmation of SARS-CoV-2 infection.

[Procedural control]

The internal procedural control zone is in the end of the membrane of the test device. STANDARD F Analyzers read the fluorescence signal of the internal procedural control zone and decide whether the result is valid or invalid.

2. The invalid result denotes that the fluorescence signal is not within the preset range. If the screen of STANDARD F analyzers shows 'Invalid Device', turn off and turn on the analyzer again and re-test with a new test device.

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BIBLIOGRAPHY

- Clinical management of severe acute respiratory infection when novel coronavirus(nCoV) infection is suspected. Interim guidance. WHO.2020
- Diagnostic detection of Wuhan coronavirus 2019 by real-time RT-PCR. 2020
- Diagnosis and treatment of pneumonia caused by new coronavirus (Novel coronavirus) National Health Commission. 2020
- Guo L et al. Profiling Early Human Response to Diagnose Novel Coronavirus Disease (COVID-19). Clinical Infectious Disease. 2020
- Zhou J et al. Antibody responses to SARS-CoV-2 in patients of novel coronavirus disease 2019. Clinical Infectious Disease. 2020

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