

STANDARD F
Dengue IgM/IgG FIA

STANDARD™ F Dengue IgM/IgG FIA
PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST SD BIOSENSOR

KIT CONTENTS



MATERIALS REQUIRED BUT NOT PROVIDED

- STANDARD F Analyzer
- Timer

SPECIMEN COLLECTION AND PREPARATION

■ Serum

- Collect the whole blood into the commercially available plain tube, NOT containing anti-coagulants such as heparin, EDTA or sodium citrate by venipuncture and leave to settle for 30 minutes for blood coagulation and then centrifuge blood to get serum specimen of supernatant.
- If serum in the plain tube is stored in a refrigerator at 2 ~ 8°C / 36 ~ 40°F, the specimen can be used for testing within 1 week after collection. Using the specimen in the long-term keeping more than 1 week can cause non-specific reaction. For prolonged storage, it should be at below -40°C / -40°F.
- It should be brought to room temperature prior to use.

■ Plasma

- Collect the venous blood into the commercially available anti-coagulant tube such as heparin, EDTA or sodium citrate by venipuncture and centrifuge blood to get plasma specimen of supernatant.
- If plasma in an anti-coagulant tube is stored in a refrigerator at 2 ~ 8°C / 36 ~ 40°F, the specimen can be used for testing within 1 week after collection. Using the specimen in the long-term keeping more than 1 week can cause non-specific reaction. For prolonged storage, it should be at below -40°C / -40°F.
- It should be brought to room temperature prior to use.

■ Whole blood

■ Capillary whole blood

- Capillary whole blood should be collected aseptically by fingertip.
- Clean the area to be lanced with an alcohol swab.
- Squeeze the end of the fingertip and pierce with a sterile lancet.
- Collect the capillary whole blood to the black line of the disposable dropper for the testing.
- The capillary whole blood must be tested immediately after collection.

■ Venous Whole blood


- Collect the venous whole blood into the commercially available anti-coagulant tube such as heparin, EDTA or sodium citrate by venipuncture.
- If venous whole blood in an anti-coagulant tube is stored in a refrigerator at 2 ~ 8°C / 36 ~ 40°F, the specimen can be used for testing within 1 ~ 2 days after collection.
- Do not use hemolyzed blood specimens.

- Anticoagulants such as heparin, EDTA or sodium citrate do not affect the test result.
- As known relevant interference, hemolytic specimen, rheumatoid factors contained specimen and lipemic, icteric specimen can lead to impact the test results.
- Use separate disposable materials for each specimen in order to avoid cross-contamination which can cause erroneous results.

TEST PROCEDURE

■ Preparation

- Allow kit components and collected specimen to room temperature at least 30 minutes before starting the test.
- Carefully read instructions for using the STANDARD F Dengue IgM/IgG FIA.
- Check the expiry date at the back of the foil pouch. Use another lot, if expiry date has passed.
- Open the foil pouch, and check the test device and the desiccant in the foil pouch.



■ Analysis of specimen

STANDARD TEST™ mode

- Prepare a STANDARD F Analyzer. Take the test device out of the foil pouch and place it on a flat and dry surface. Write patient information on the label of test device. Select the 'Standard Test' mode according to the analyzer's manual as below.

STANDARD F2400 Analyzer	Worksheet → "Run Test" → Scan or type patient ID and/or operator ID
STANDARD F200 Analyzer	Standard Test mode → Insert patient ID and / or operator ID on the analyzer

- Insert the test device to the test slot of the analyzer. When inserting the test device to the analyzer, the analyzer will read the barcode data, and check the test device is valid.



- Collect the 10 µl of serum/plasma/whole blood to the black line of a the STANDARD Exti Tube.



- Add the collected serum/plasma/whole blood to the specimen well of the test device.



- Add 3 drops of assay diluent to the assay diluent well of the test device.



PERFORMANCE CHARACTERISTICS

■ CLINICAL PERFORMANCE

43 of Dengue IgM/IgG positive confi rmed specimens and 384 of Dengue IgM/IgG negative confi rmed specimens were tested internally. For STANDARD F Dengue IgM/IgG FIA, the result shows the 98.0%/42/43 sensitivity and 383/384 (99%) specificity data for plasma specimens.

Clinical Sensitivity	Reference test (ELISA)			
	Dengue IgM/IgG	POS	NEG	TOTAL
STANDARD™ F Dengue IgM/IgG FIA	POS	42	0	42
	NEG	1	0	1
	TOTAL	43	0	43
Sensitivity: 42/43 (98%)				

Clinical Specificity	Reference test (ELISA)			
	Dengue IgM/IgG	POS	NEG	TOTAL
STANDARD™ F Dengue IgM/IgG Test	POS	0	1	1
	NEG	0	183	183
	TOTAL	0	184	184
Specificity: 183/184 (99%)				

■ ANALYTICAL PERFORMANCE

1. Analytical Sensitivity (Limit of Detection)

The analytical sensitivity (Limit of Detection [LOD]) of STANDARD™ F Dengue IgM/IgG FIA are shown below.

STANDARD™ F Dengue IgM/IgG FIA	
Antibody	Limit of Detection(ng/mL)
Dengue IgM	2.0
Dengue IgG	3.9

- After applying the specimen, immediately press the "TEST START" button.



- The analyzer will automatically display the test result after 15 minutes.



READ ONLY™ mode

STANDARD F200 Analyzer

- Take the test device out of the foil pouch and place it on a flat and dry surface. Write a specimen information on the label of test device.
- Collect the 10 µl of serum/plasma/whole blood to the black line of a the STANDARD Exti Tube.



- Add the collected serum / plasma / whole blood to the specimen well of the test device.



- Add 3 drops of assay diluent into the assay diluent well of the test device.



- Incubate the test device for 15 minutes outside of the analyzer. Incubation must not be more than 20 minutes.



- Prepare the STANDARD F Analyzer and set the "READ ONLY" mode following the instructions in the manual.
- Insert the test device to the test slot of the analyzer. When inserting the test device to the analyzer, the analyzer will automatically scan and display the test results.



INTERPRETATION OF TEST RESULTS

Result	COI (Cut-off index) value	Interpretation
Positive	COI ≥ 1.0	Positive for Dengue IgM/IgG antibody
Negative	COI < 1.0	Negative for Dengue IgM/IgG antibody
Invalid	Not show the COI value	Retest should be performed

Results should be considered in conjunction with the clinical history and other data available to the physician.

The Analyzer's test result of a specimen is given either as Positive (+) / Pos(+) or Negative(-) / Neg(-) with a COI (cut-off index) value. Cut-off index (COI) is based on the ratio of assay signal to cut-off value.

EXPLANATION AND SUMMARY

■ Introduction

Dengue viruses, transmitted by Aedes aegypti and Aedes albopictus mosquitoes, are widely distributed throughout the tropical and subtropical areas of the world. There are four known distinct serotypes of dengue virus (DENV-1, DENV-2, DENV-3 and DENV-4). Rapid and reliable tests for primary and secondary infections of Dengue are essential for patient management. An infected person experiences the acute symptoms of Dengue when there is a high level of the virus in the bloodstream. As the immune response fights the Dengue infection, the person's cells begin producing IgM and IgG antibodies that are released in the blood and lymph fluids, where they recognize and neutralize the Dengue virus and viral molecules such as the Dengue non-structural protein 1 (NS1) antigen.

■ Intended use

STANDARD F Dengue IgM/IgG FIA is a fluorescence immunoassay for the detection of IgM/IgG antibodies against Dengue virus in human serum, plasma, or whole blood specimens. This test kit is for in vitro use only. This is intended for professional use only for an initial screening test. Test results of this kit have to analyze with appropriate analyzer, STANDARD F Analyzer, manufactured by SD BIOSENSOR.

■ Test principle

STANDARD F Dengue IgM/IgG FIA has "M", "G" test lines and "C" control line. Monoclonal anti-human IgM and monoclonal anti-human IgG are immobilized at two individual test lines respectively (M, G line) on the nitrocellulose membrane. Inactivated Dengue virus in the antigen pad and europium conjugated antibodies (monoclonal anti-Dengue Env-Ep and monoclonal anti-Dengue NS1-Ep) in the conjugation pad release by adding assay diluents and react with anti-Dengue IgM or IgG in patient specimen. If human anti-Dengue IgM or IgG exist in patient serum, complexes with anti-human IgM/IgG on the test lines, human IgM/IgG in patient specimen, inactivated Dengue virus in the antigen pad, and europium conjugated antibodies in the conjugation pad make fluorescence signal. The intensity of the fluorescence light generated on the membrane is scanned by the STANDARD F Analyzer manufactured by SD BIOSENSOR. STANDARD F Analyzer can analyze the presence of the analyte in the clinical specimen by processing the results using pre-programmed algorithms and display the test result on the screen.

KIT STORAGE AND STABILITY

Store the kit at 2 ~ 30°C / 36 ~ 86°F, out of direct sunlight. Kit materials are stable until expiration date printed on the outer box. Do not freeze the kit.

WARNINGS AND PRECAUTIONS

- Do not re-use the test kit.
- Use the STANDARD F Dengue IgM/IgG FIA at 15 ~ 32°C / 59 ~ 90°F and 10 ~ 90% RH.
- Do not use the kit if the pouch is damaged or the seal is broken.
- Do not smoke, drink or eat while handling specimen.
- Wear personal protective equipment, such as gloves and lab coats when handling kit reagents. Wash hands thoroughly after you experiment.
- Clean up spills thoroughly using an appropriate disinfectant.
- Handle all specimens as if they contain infectious agents.
- Observe established precautions against microbiological hazards throughout testing procedures.
- Dispose of all specimens and materials used to perform the test as bio-hazard waste. Laboratory chemical and biohazard wastes must be handled and discarded in accordance with all local, state, and national regulations.
- The barcode of the test device is used by analyzer to identify the type of test being run and to identify the individual test device so as to prevent to a second read of the test device by the same analyzer.
- Immediately use the test device after taking out of a foil pouch.
- As the detection reagent is a fluorescent compound, no visible results will form on the test device. The STANDARD F Analyzers authorized by SD BIOSENSOR must be used for result interpretation.
- Improper specimen collection, handling or transport may yield inaccurate results.
- Do not write on the barcode or damage the barcode of the test device.

LIMITATION OF TEST

- The contents of this kit are to be used the qualitative detection of anti-Dengue IgM/IgG from blood specimens of the symptomatic patients.
- Failure to follow the test procedure or improper specimen collection may adversely affect test performance or invalidate the test result.
- For more accuracy of immune status, additional follow-up testing using other laboratory methods is recommended.
- A negative test result may occur if the level of antibody in a specimen is below the detection limit of the test or if the specimen was collected, transported, or stored improperly.
- Negative test results do not rule out possible other infections.
- Positive test results do not rule out coinfection with other pathogens.

QUALITY CONTROL

■ Internal procedural control

- The internal procedural control zone is on the membrane of the test device. STANDARD F Analyzers read the fluorescence signal of the internal procedural control zone and decide whether the result is valid or invalid. The invalid result is due to the test that the fluorescence signal is not within the pre-set range. If the screen displays STANDARD F Analyzers shows 'Invalid Device', turn off and turn on the analyzer again and re-test with the new test device.

No.	Type	Potential Cross-reacting Microorganism/Virus	Concentration
1		Adinetobacter baumannii	
2		Bacteroides fragilis	
3		Bordetella pertussis	
4		Candida albicans	
5		Chlamydia pneumoniae	
6		Escherichia coli	
7	Bacteria	Haemophilus influenzae	2.0x10 ⁷ cfu/mL
8		Haemophilus influenzae	
9		Kingella kingae	
10		Listeria monocytogenes	
11		Lactobacillus plantarum	
12		Legionella pneumophila	
13		Moraxella catarrhalis	
14		Mycobacterium avium	
15		Mycobacterium tuberculosis	
16		Mycoplasma pneumoniae	
17		Neisseria gonorrhoeae	
18		Neisseria meningitidis	
19		Neisseria mucosa	
20		Neisseria sicca	
21		Peptostreptococcus anaerobius	
22		Prevotella oralis	
23		Propionibacterium acnes	
24	Bacteria	Proteus mirabilis	2.0x10 ⁷ cfu/mL
25		Pseudomonas aeruginosa	
26		Serratia marcescens	
27		Staphylococcus aureus	
28		Staphylococcus epidermidis	
29		Streptococcus mutans	
30		Streptococcus pneumoniae	
31		Streptococcus pyogenes	
32		Streptococcus salivarius	
33		Streptococcus sanguis	
34		Vibrio parvula	

35		Adenovirus 3	
36		Adenovirus 4	
37		Adenovirus 5	
38		Adenovirus 11	
39		Coronavirus 229E	
40		Coronavirus OC43	
41		Cytomegalovirus AD-169	
42		Cytomegalovirus Towne	
43		Echovirus Type 3	
44		Enterovirus	
45	Viruses	Herpes Simplex virus 1	2.0x10 ⁷ TCID ₅₀ /mL
46		Herpes Simplex virus 2	
47		HSV Type 1	
48		Human Coronavirus OC43	
49		Human Metapneumovirus A1	
50		Human Metapneumovirus A2	
51		Human Metapneumovirus B1	
52		Human Metapneumovirus B2	
53		Human Parainfluenza	
54		Influenza A H1N1 (Denver/1/57)	
55		Influenza A H1N1 (FM/1/47)	

56	Viruses	Influenza A H1N1 (Mexico/4108/2009)	2.0x10 ⁷ TCID ₅₀ /mL
57		Influenza A H1N1 (New Jersey/8/76)	
58		Influenza A H1N1 (PR/8/34)	
59		Influenza A H3N2	
60		Influenza B Hong Kong	
61		Influenza B Panama	
62		Influenza C Taylor	
63		Measles virus	
64		Mumps virus	
65		Parainfluenza virus 1	
66		Parainfluenza virus 2	
67		Parainfluenza virus 3	
68		Parainfluenza virus 4A	
69		Parainfluenza virus 4B	
70		Rhinovirus Type 2	

No.	Interfering Substances	Concentration
1	Acetaminophen	20 mg/mL
2	Acetylsalicylic acid	20 mg/mL
3	Albuterol	20 mg/mL
4	Betamethasone	10 µg/mL
5	Bilirubin	200 µg/mL
6	Budesonide	10 µg/mL
7	Chlorpheniramine	5 mg/mL
8	Dextromethorphan	10 mg/mL
9	Diphenhydramine	4 mg/mL
10	Fluticasone	500 ng/mL
11	Fluticasone	500 ng/mL
12	Guaiacol	30 mg/mL
13	Hemoglobin	10 mg/mL
14	Homeopathic Allergy MedDore	20 mg/mL
15	Ibuprofen	20 mg/mL
16	Maclin	10 mg/mL
17	Oxymetazoline	0.05 mg/mL
18	Phenylephrine	10 mg/mL
19	Ribavirin	1 µg/mL
20	Rimantadine	500 µg/mL
21	Synagis	4 µg/mL
22	Tobramycin	500 ng/mL

3. Interference Substances

The following substances do not interfere the test result up to the indicated concentrations;

23	Triamcinolone	500 ng/mL
24	Triglycride	9 mg/mL

BIBLIOGRAPHY

1. Dengue guidelines for diagnosis, treatment, prevention and control, World Health Organization, New Edition 2009.
2. Kilis SC, Nimmanitya S, Nisalak A, Burke DS, Evidence that maternal dengue antibodies are important in the development of dengue hemorrhagic fever in infants, Am J Trop Med Hyg Jan, 38(2):411-419, 1988.
3. Dengue haemorrhagic fever: Diagnosis, treatment, prevention and control, World Health Organization 2nd Edition, 1997.
4. Luthfi D. et al., Serological differentiation of infections with dengue virus serotypes 1 to 4 by using recombinant antigens, J Clin Microbiol, 40(11):4317-4320, 2002.
5. Matthew T. R. et al. Dengue virus pirates human platelets, Blood, 126(3):286-287, 2015.
6. Guzman M. G. et al. Dengue: A continuing global threat, Nat Rev Microbiol, 8:57 -516, 2010.

SYMBOL											
	Reference number		Caution		Use by		Batch code		Consult Instructions for Use		Do not re-use
	In vitro Diagnostics		Note		Manufacturer		Date of manufacture		Contains Sufficient for <= 5 Tests		Keep away from sunlight
	Indicate that you should keep the product dry		To indicate the temperature limitations in which the transport package has to be kept and handled.				This product fulfills the requirements of the European Directive 98/79/EC.			Do not use if packaging is damaged	

EC | REP

Authorized Representative
MT Promedk Consulting GmbH
Ermst-Heckel-Straße 7 66386 St. Ingbert Germany
Phone : +49 6894 581020, Fax : +49 6894 581021

 **Manufactured by SD Biosensor, Inc.**

Head office : C-40&5th 16, Daejeong-daero 3556beon-gil, Yeongtong-gu, Suwon-si, Gyeonggi-do, 16690, REPUBLIC OF KOREA
Manufacturing site : 74, Osongsaengmyeong 4-ro, Osong-eup, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, 28161, REPUBLIC OF KOREA

Please contact us for any complaints/inquiries/suggestions via email (sales@sdbiosensor.com), phone (+82-31-300-0400) or website (www.sdbiosensor.com).