STANDARD F Tnl Pro FIA

STANDARD™ F Tnl Pro FIA

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

KIT CONTENTS

MATERIALS REQUIRED BUT NOT PROVIDED

- STANDARD F Analyzer

SPECIMEN COLLECTION AND PREPARATION

Whole blood

Venous whole blood

- Collect the venous whole blood into the commercially available EDTA tube by venipuncture.
- It is recommended that collected venous whole blood specimens are used immediately. If venous whole blood in an anticoagulant tube is stored at room temperature (15-25 °C / 59-77 °F) and refrigerator conditions (2-8 °C / 36-46 °F), the specimen can be used for testing within 8 hours after collection.
- Do not use hemolyzed blood specimens.


Serum

- Collect the whole blood into the commercially available plain tube. NOT containing anti-coagulants such as heparin, EDTA or sodium citrate by venipuncture and leave to settle for 30 minutes for blood coagulation and then centrifuge blood to get serum specimen of supernatant.
- Serum specimens may be stored at room temperature (15-25 °C / 59-77 °F) and refrigerator conditions (2-8 °C / 36-46 °F) for 1 day prior to testing.

TEST PROCEDURE


Preparation

- Allow kit components and collected specimen to 15-30 °C / 59-86 °F at least 30 minutes before starting the test.
- Carefully read instructions for the STANDARD F Tnl Pro FIA.
- Check the expiry date at the back of the foil pouch. Use another lot if expiry date has passed.
- Open the foil pouch and check the test device in the foil pouch.

Analysis of control

'STANDARD TEST' mode	STANDARD F200 and F2400 Analyzer
1. Take the test device out of the foil pouch and place it on a flat and dry surface. Write patient information on the label of test device.	Workplace' → 'Run Test' → Insert patient ID and / or operator ID on the analyzer
2. Insert the test device to the test slot of the analyzer. When inserting the test device to the analyzer, the analyzer read the barcode data, and check the test device is valid.	
3. Using a disposable dropper, collect the 100 µl of the specimen (whole blood or serum) to the black line of the disposable dropper. Add the collected specimen into the extraction buffer and mix thoroughly.	
4. Apply the 100 µl of specimen mixture to the specimen well of the test device.	'Standard Test' mode → Insert patient ID and / or operator ID
5. After applying the control mixture, immediately press the "TEST START" button.	
6. The analyzer will automatically display the test result after 10 minutes.	

PERFORMANCE CHARACTERISTICS

Analytical performance

1) Accuracy (Method comparison)
Results of the serum / whole blood measured by STANDARD™ F Tnl Pro FIA meet the acceptance criteria as below;

Regression Analysis		
	SERUM	WHOLE BLOOD
Slope	0.9704	1.0022
Y-intercept	47.971	7.425
R	0.9955	0.9945
R ²	0.9910	0.9890
n	210	210

INTERPRETATION OF TEST RESULTS

The reportable range is 10 ng/L to 20,000 ng/L. If result is below 10 ng/L, it will be reported as “ $\leq$ 10 ng/L”. If result is above 20,000 ng/L, it will be reported as “ $\geq$ 20,000 ng/L”.

Measuring range	0.01 - 20 ng/mL
Cut off value for IMI	< 0.05 ng/mL

The default results are measured in ng/L. Should an alternate result unit be desired, the user can change unit through the “Setting” menu of the analyzer or manually calculate using the following conversion factors:

(Concentration in ng/L) x (0.001) = ng/mL = µg/L
(Concentration in ng/L) x (1) = pg/mL

EXPLANATION AND SUMMARY

Introduction

The Troponin (TnI, one of the subunits of the troponin regulatory complex, binds to actin and inhibits interactions between actin and myosin. Troponin complex plays an important role in the regulation of skeletal and cardiac muscle contraction. Troponin complex consists of three subunits: Troponin T (TnT), Troponin I (TnI) and Troponin C (TnC). The subunits are held together by non-covalent interactions. TnI is the troponomyosin-binding subunit that regulates interaction between the troponin complex and the thin filament. The TnI subunit is responsible for inhibiting actomyosin activation at low intracellular Ca²⁺ concentrations. The third subunit, TnC, binds Ca²⁺ ions during the excitation of the muscle and changes the conformation of troponin complex, thus enabling the formation of actomyosin complex and the subsequent muscle contraction. Cardiac Troponin I (cTnI) has been known as reliable markers of cardiac cell death. It is now widely used for the diagnosis of acute myocardial infarction (AMI), postsurgery myocardium trauma, chemotherapy cardiotoxicity, as well as some other diseases related to cardiac muscle injury. Being expressed cTnI only in cardiac tissue, cardiac TnI have been the preferred biomarkers for myocardial infarction for a long time. STANDARD F Tnl Pro FIA, employing immunofluorescent detection system with STANDARD F analyzer, provides significantly fast, easy and accurate system to identify cardiac Troponin I (cTnI) in serum and whole blood as an aid in the diagnosis of myocardial infarction and enables supportive treatment decisions.

Intended use

STANDARD F Tnl Pro FIA is a fluorescence immunoassay for the quantitative determination of cardiac Troponin I (cTnI) levels in human serum and whole blood using STANDARD F Analyzers, manufactured by SD BIOSENSOR. This test is for *in vitro* diagnostic use and is intendeded for use as an aid in the screening and monitoring of acute myocardic infarction (MI).

Test principle

STANDARD F Tnl Pro FIA uses antibodies that specifically bind to the TnI protein. When a human specimen is applied into the specimen well of the test device, the specimen migrates through the membrane from the specimen well. If Troponin I (TnI) is present, it will be bound to the detector coupled to europium microparticle that migrates through the membrane. The europium microparticle complex will be captured by capture antibody on the test line where it is detected by STANDARD F Analyzer. If it is not present, the europium microparticle will not be trapped by the capture antibody nor detected by STANDARD F Analyzer. The intensity of the fluorescence light generated on the membrane is scanned by the STANDARD F Analyzer. STANDARD F Analyzer can analyze the presence of the TnI in the clinical specimen by processing the results using pre-programmed algorithms and display the test result on the screen.

KIT STORAGE AND STABILITY

Store the kit at 2-30 °C / 36-86 °F out of direct sunlight. Kit materials are stable until expiration date printed on the outer box. Do not freeze the kit.

WARNINGS AND PRECAUTIONS

- STANDARD F Tnl Pro FIA is for *in vitro* diagnostic use only.
- Carefully follow the instructions and procedures described in these instructions for use before testing.
- STANDARD F Tnl Pro FIA should be used with STANDARD F Analyzer.
- STANDARD F Tnl Pro FIA should remain in its original sealed pouch until ready to use. Do not use the test device if pouch is damaged or the seal is broken.
- STANDARD F Tnl Pro FIA test device is single use only. Do not re-use it.
- Do not use hemolyzed specimens or frozen specimens.
- Do not use any artificial materials.
- Place the analyzer on a flat surface when in use.
- Wash your hands in warm and soapy water. Rinse well and dry completely before testing.
- Discard the used test kit according to the proper method.
- Desiccant in foil pouch is to absorb moisture and keep humidity from affecting products. If the moisture indicating Desiccant beads change from yellow to green, the test device in the pouch should be discarded.
- Check the expiration date printed at the pouch or package.
- Check the volume (100 µl) of extraction buffer.
- Use the STANDARD F Tnl Pro FIA at 15-30 °C / 59-86 °F.
- All kit components must be at room temperature 30 minutes before running the assay.
- Do not write on the barcode or damage the barcode of the test device.

LIMITATION OF TEST

- The test should be used for the detection of Troponin I levels in human whole blood and serum specimens.
- Failure to follow the test procedure and interpretation of test results may adversely affect test performance and/ or produce invalid results.
- Invalid results may occur if a poor-quality specimen is obtained.
- The test result must always be evaluated with other clinical data available to the physician.

QUALITY CONTROL

Calibration

The calibration set test of STANDARD F Analyzer should be conducted according to the analyzer's manual.


[When to use calibration set]

- Before using the analyzer for the first time.
- When you drop the analyzer.
- Whenever you do not agree with the final result.
- When you want to check the performance of an analyzer and test device.

[How to use calibration set]

Calibration set test is a required function that ensures optimal performance by checking the internal analyzer optics and functions.

- Select the 'Calibration' in main menu.
- The specific calibration set is included with the analyzer.
- Insert the CAL-1 for white calibration, CAL-2 for UV LED calibration, and CAL-3 for RGB LED calibration in sequence.

Internal quality control

- The internal procedural control zone is on the membrane of the test device. STANDARD F Analyzers read the fluorescence signal of the internal procedural control zone and decide whether the result is valid or invalid.
- The invalid result denotes that the fluorescence signal is not within the pre-set range. If the screen of STANDARD F Analyzer shows 'Invalid', turn off and turn on the analyzer again and re-test with a new test device.

External quality control

Quality control testing should be run to check the performance of STANDARD F Tnl Pro FIA and STANDARD F Analyzer. STANDARD F Tnl Control (manufactured by SD BIOSENSOR can be used for the external quality control. Control test should be conducted in accordance with the instructions of STANDARD F Tnl Control.

External quality control test should be run:

- once for new lot
- once for each untrained operator
- as required by test procedures in the Instructions for use of STANDARD F Tnl Control and in accordance with local, state and federal regulations or accreditation requirements

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


3) Precision
Precision of STANDARD™ F Tnl Pro FIA has been evaluated through the within-run and day-to-day study. All results met the acceptance criteria (CV≤15%).

Within-run Study										
Specimen type	TnI (ng/L)	Level 1			Level 2			Level 3		
		Site 1	Site 2	Site 3	Site 1	Site 2	Site 3	Site 1	Site 2	Site 3
Serum	Ref. value	n	100	100	100	100	100	100	100	100
	Ref. value	70	70	70	6500	6210	6260	13420	13750	13540
	Mean	62.13	62.36	62.62	6476	6447	6478	14403	14288	14503
	CV	7.4%	7.2%	7.5%	8.6%	8.5%	8.9%	8.1%	8.5%	7.51%
	n	100	100	100	100	100	100	100	100	100
	Ref. value	70	70	70	6360	6580	6390	14100	14130	14230
Whole blood	Ref. value	62.46	62.68	62.31	6532	6501	6590	14441	14249	14278
	Mean	62.46	62.68	62.31	6532	6501	6590	14441	14249	14278
	CV	8.0%	7.4%	7.4%	9.0%	9.7%	8.5%	8.21%	8.41%	7.71%

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