STANDARD **M10****Flu/RSV/SARS-CoV-2 Fast**

STANDARD™ M10 Flu/RSV/SARS-CoV-2 Fast
PLEASE READ INSTRUCTIONS CAREFULLY BEFORE
YOU PERFORM THE TEST

**Intended Purpose**

STANDARD M10 Flu/RSV/SARS-CoV-2 Fast test, performed on the STANDARD M10 system, is an automated multiplexed real-time RT-PCR assay intended for the simultaneous *in vitro* qualitative detection and differentiation of severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2), influenza A, influenza B, and/or respiratory syncytial virus (RSV) RNA in nasopharyngeal swab specimens collected from individuals presenting signs and symptoms of respiratory tract infection.

The test is intended to aid in the diagnosis of COVID-19, influenza, and/or RSV infections and it is for professional use only and is intended to be used by clinical laboratory professionals or healthcare professionals who have received general professional training appropriate to their role. The test is designed for use in clinical laboratories or near-patient testing environments, such as healthcare facilities and other adequately equipped healthcare settings.

Positive results are indicative of the presence of Influenza A, Influenza B, RSV and/or SARS-CoV-2 RNA. Clinical correlation with patient history and other diagnostic information is necessary to determine patient infection status.

Negative results must be combined with clinical observations, patient history, and epidemiological information.

Summary and Explanation

Acute respiratory infection can be caused by a variety of viruses and bacteria, including recently introduced SARS-CoV-2. Acute respiratory infection of SARS-CoV-2 outbreak in Wuhan, China has widely spread out into the world since 2019. Common signs of a person infected with SARS-CoV-2 include respiratory symptoms, fever, cough, shortness of breath, and dyspnea. In more severe cases, the infection generates pneumonia, acute respiratory syndrome, kidney failure, or even death.

Influenza, commonly known as the flu, is a highly contagious and acute viral infection of the respiratory tract caused by an influenza virus. Three types of influenza viruses affect people, called Type A, Type B, and Type C. Type A viruses are the most prevalent and are associated with most serious epidemics. The clinical symptoms

caused by influenza A virus infection are more severe than those caused by influenza B virus. Both type A and B viruses can circulate simultaneously, but usually one type is dominant during a given season. Influenza can take on a variety of appearances, ranging from isolated respiratory findings that resemble the common cold, to severe pneumonia requiring hospitalization and death.

RSV (Respiratory syncytial virus) is an enveloped, negative-sense RNA virus belonging to the Paramyxoviridae family.

It is prevalent worldwide and tends to cause seasonal outbreaks, particularly during winter months. The virus infects the cells lining the respiratory tract, leading to inflammation and the production of excess mucus.

RSV infection can result in a range of respiratory symptoms, from mild cold-like illness to more severe conditions such as bronchiolitis or pneumonia.

Due to its potential for causing serious illness and contributing to significant healthcare burden globally, the rapid and accessible detection of RSV is essential for timely diagnosis, infection control, and patient management.

Respiratory infections caused by SARS-CoV-2, influenza, and RSV have similar transmission routes and symptoms. Therefore, the use of assays capable of simultaneous detection and differentiation enables efficient control of viruses, appropriate treatment, and prevention of widespread outbreaks.

■ Cartridge Description

STANDARD M10 Flu/RSV/SARS-CoV-2 Fast is a molecular *in vitro* diagnostic assay that aids in the simultaneous detection and differentiation of Influenza A, Influenza B, RSV and SARS-CoV-2 virus RNA based on nucleic acid amplification technology, RT-PCR. STANDARD M10 Flu/RSV/SARS-CoV-2 Fast cartridge contains viral RNA extraction buffers and RT-PCR reagents for the *in vitro* qualitative detection of Influenza A, Influenza B, RSV and SARS-CoV-2 virus RNA in nasopharyngeal swab specimens.

Figure 1. Layout of STANDARD M10 Flu/RSV/SARS-CoV-2 Fast cartridge

Principle of the Procedure

STANDARD M10 Flu/RSV/SARS-CoV-2 Fast test is an automated *in vitro* diagnostic test for the qualitative detection of nucleic acid from Flu A, Flu B, RSV and SARS-CoV-2. STANDARD M10 Flu/RSV/SARS-CoV-2 Fast test is performed on STANDARD M10 system.

STANDARD M10 system automates and integrates sample preparation, nucleic acid extraction, reverse transcription polymerase chain reaction (RT-PCR), and detection of the target sequences in various specimens using molecular diagnostic assays. The system consists of the STANDARD M10 Module and the STANDARD M10 Console with preloaded software for running tests and viewing the results. The system requires the use of single-use disposable cartridges that hold the RT-PCR reagents and host the RT-PCR process. Because the cartridges are self-contained, cross-contamination between samples is minimized. For a full description of the systems, see the STANDARD M10 system User Manual.

STANDARD M10 Flu/RSV/SARS-CoV-2 Fast test includes reagents for the detection of RNA from Flu A, Flu B, RSV and SARS-CoV-2 in nasopharyngeal swab specimens. The cartridge is present to control for adequate processing of the sample, nucleic acid extraction and RT-PCR. The sample is first lysed under highly denaturing conditions to inactivate RNase and to ensure isolation of intact viral RNA. Buffering conditions are then adjusted to provide optimal binding of the RNA to the glass fibers. The RNA binds to the glass fibers, and contaminants are efficiently washed away several times using a wash buffer. High-quality RNA is eluted with elution buffer, ready for amplification. The purified nucleic acids are mixed with the lyophilized PCR master mix beads. And it is denatured under the customized temperature conditions to carry out effective multiplex real-time RT-PCR and amplified according to the PCR conditions.

The table below indicates which target is designed to be detected by three channels: FAM, HEX, Cy5.

Table 1. Fluorescent channel of each target Pathogen

Pathogen	Target gene	Channel
Influenza	Flu A (M gene, PB2 gene)	FAM
	Flu B (NS1 gene, M gene)	HEX
	Internal control (IC)	Cy5
RSV	RSV (M gene, N gene)	HEX
	Internal control (IC)	Cy5
SARS-CoV-2	ORF1ab gene, E gene	FAM
	N gene	HEX
	Internal control (IC)	Cy5

Active Ingredient

- Guanidine Thiocyanate 0.47g
- Taq DNA Polymerase 646U
- Reverse Transcriptase 1149U
- Target primers and probes 2.5 - 30pmole

Materials Provided

STANDARD M10 Flu/RSV/SARS-CoV-2 Fast kit contains sufficient reagents to process 10 specimens or external quality control samples.

Table 2. Contents of STANDARD M10 Flu/RSV/SARS-CoV-2 Fast kit

	Contents	Quantity	Usage in each reaction
1	Cartridge	10	1ea
2	STANDARD Fixed volume dropper (300 µL)	11	1ea
3	Instructions for Use	1	-

Storage and Handling

Store STANDARD M10 Flu/RSV/SARS-CoV-2 Fast kit at 2-28 °C (36-82 °F). If the cartridge has been refrigerated, perform the test after stabilizing it for 30minutes at 20-28 °C (68-82 °F).

Do not store the cartridge upside-down. If stored upside-down, it remains stable for up to 14 days. Do not remove the Safety Clip of the cartridge and do not press the cartridge until actual use. Do not use a cartridge that has leaked or is wet. Under these conditions, cartridges can be stored until the expiration date printed on the packaging.

Materials Required but Not Provided

- STANDARD™ M10 system with User Manual
 - At least one STANDARD M10 Console (Cat. No. 11M1011) and one STANDARD M10 Module (Cat. No. 11M1012)
 - * For STANDARD M10 Console S/W is V 001.010 or higher version
 - STANDARD M10 Calibration Kit (Cat. No. 11CAL10H)
- Sample collection tools
 - Copan eNAT (6U074S01, 606CS01P)
 - Copan UTM (305C)
 - Noble bio CTM (UTNFS-3B-2)
- Sample transfer pipettes (Optional)
 - Micropipette with filter tips
- Biohazard container
- PPE (Personal Protective Equipment)
- External Quality Control Material
 - AccuPlex™ SARS-CoV-2, Flu A/B and RSV Molecular Controls Kit (SeraCare, 0505-0260)

Warnings and Precautions

1. This kit is for *in vitro* diagnosis.
2. Please read the Instructions for Use carefully before starting the test.
3. Improper specimen collection, transfer, storage, and processing may cause erroneous test results.
4. Do not remove the Safety Clip of the cartridge before use.
5. Do not press the cartridge until actual use.
6. Press (punch) the cartridge and ensure that the cartridge lid (Top cover) is fully engaged with the cartridge body before dispensing the sample.
7. Insert the tip or the supplied dropper deeply into the sample port to ensure that the sample is dispensed accurately.
8. Do not use a cartridge that has leaked or is wet.
9. Do not use the kit after its expiration date.
10. Do not use the product if more than 3 months have passed since opening the foil pouch. The specified stability after opening only applies when the cartridge is stored at 2–28°C.
11. Do not shake, tilt, or invert the cartridge especially after pressing the cartridge to punch the seal. It may yield invalid or false test results.
12. Do not use a cartridge with a damaged barcode label.
13. Do not reuse processed cartridges.
14. Biological specimens, transfer devices, and used cartridges should be considered capable of transmitting infectious agents requiring standard precautions.
15. This test involves the extraction of viral RNA and PCR amplification; therefore, precautions must be taken to avoid contamination in both clinical laboratories and near-patient testing (NPT) environments equipped with appropriate medical infrastructure and operated by trained healthcare professionals. Regular environmental monitoring and adherence to contamination control procedures are recommended. In case of spills or exposure, decontaminate the area immediately using an effective disinfectant such as bleach or 70% alcohol. If bleach is used, follow with sterile water to prevent damage to the testing device or workspace.
16. When using this kit, procedures must be followed strictly in accordance with the instructions for use (IFU), and testing should comply with the technical requirements applicable to clinical laboratories and near-patient testing (NPT) environments.
17. Dispose of all specimens and materials used to perform the test as biohazard waste. Laboratory chemical and biohazard wastes must be handled and discarded in accordance with all local, state, and national regulations.
18. This kit is intended to be used only by trained laboratory professional and healthcare professionals, such as physicians and nurses, who have received general professional training appropriate to their role. Use by untrained or non-qualified individuals may lead to inaccurate results.
19. The test should be operated in an environment with a temperature range of 15°C to 35°C (59°F to 95°F). Operating the test outside this range may affect test accuracy and reliability.

20. Ensure the device is placed on a stable surface, away from direct sunlight, heat sources, and excessive humidity.
21. Do not use the device outside of healthcare facilities or equivalent settings.
22. Use of other than the indicated sample collection tools, not specified for use with the STANDARD M10 Flu/RSV/SARS-CoV-2 Fast, may lead to inaccurate results and is not recommended.
23. Any serious incident that has occurred in relation to this kit shall be reported to the manufacturer or its authorized representative, and to the competent authority of the Member State in which the user and/or the patient is established, in accordance with the requirements of Regulation (EU) 2017/746.

Chemical Hazards

This kit contains components classified as follows in accordance with the regulation (EC) No. 1272/2008:

* **Hazard pictogram:** 

* **Signal word:** Danger

* **Hazard statements:**

-H314 Causes severe skin burns and eye damage

-H318 Causes serious eye damage

* **Precautionary statements:**

1) **Prevention**

-P260 Do not breathe dust/fume/gas/mist/vapours/spray.

-P264 Wash hands thoroughly after handling.

-P280 Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.

2) **Response**

-P301+P330+P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

-P303+P361+P353 IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.

-P304+P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.

-P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

-P310 Immediately call a POISON CENTER or doctor/physician.

-P321 Specific treatment

-P363 Wash contaminated clothing before reuse.

3) **Storage**

-P405 Store locked up

4) **Disposal**

-P501 Dispose of contents/container in accordance with local/regional/national/international regulation

Specimen Collection, Transport, and Storage

Proper sample collection, transportation, and storage are critical to the performance of the test. Improper sample collection, inappropriate sample handling and/or transportation can lead to false results.

For comprehensive guidance, users should refer to the relevant national or local biosafety standards and specimen management procedures applicable to respiratory virus testing.

■ Nasopharyngeal swab collection procedure

1. Gently and slowly insert the nasopharyngeal swab through the nostril parallel to the back of the nasopharynx.
2. Rotate it several times to obtain secretions.
3. Remove the swab from the nasopharynx and place it into the specimen collection tube, and discard the tail.
4. Cap the specimen collection tube tightly. If testing is to be performed immediately, the sample can be loaded into the cartridge immediately at this step. For detailed test method, refer to 'Loading a sample into STANDARD M10 Flu/RSV/SARS-CoV-2 Fast Cartridge'.
5. The swab specimens in the collection tube can be stored up to 12 hours at 19-25°C (66.2-77°F), 4 days at 2-8 °C (36-46 °F), 3 months at -18°C (-0.4°F) or below. Up to 3 freeze/thaw cycles are acceptable.


Figure 2. Nasopharyngeal swab collection

PROCEDURE

■ Starting the STANDARD™ M10 system

Note

For the detailed instructions, refer to the STANDARD M10 system User Manual. If you have scanned the cartridge barcode in the STANDARD M10 and the software version is not compatible, a 'Not Supported Device' error message appears. Update the software before proceeding with the test.

1. Turn on the STANDARD M10 system.
2. Check the STANDARD M10 Console and the STANDARD M10 Module is connected and functional.

Figure 3. Power connection

3. Enter the User ID and Password on the Log In screen of the STANDARD M10 Console and click the Log In button.
4. Touch the STANDARD M10 Module to run on the Home screen.
(The door of the selected STANDARD M10 Module will automatically open for cartridge loading.)


Figure 4. Log In screen

Figure 5. Home screen, Status of M10 module

5. Enter a Patient ID by scanning the barcode or using virtual keyboard on the M10 Console screen. (Patient ID is optional. You can turn off the Patient ID option from the 'Settings'.)
6. Enter a Sample ID by scanning the barcode of the specimen or using virtual keyboard on the M10 Console screen. Make sure that the specimen tube cap is firmly closed when scan the ID barcode printed on the specimen tube. (For quality control test, tick the QC check box.)

Figure 6. Entering Sample ID

Figure 7. Scanning a cartridge

- Scan STANDARD M10 Flu/RSV/SARS-CoV-2 Fast cartridge to be used. The STANDARD M10 Console automatically recognizes the assay to be run based on the cartridge barcode.

Note If you have scanned the cartridge barcode in the STANDARD M10 and the expiration date has expired, An 'Expired Device' error message appears. Check validity period and test with unexpired cartridges.

■ Loading a sample into STANDARD M10 Flu/RSV/SARS-CoV-2 Fast cartridge

Caution If the cartridge has been refrigerated, perform the test after stabilizing it for 30 minutes at 20-28 °C (68-82 °F). Once the sample has been loaded into the cartridge, start the test as soon as possible. (within 30 minutes).

Note False negative results may occur if insufficient sample is added into the cartridge.

- Remove the Safety Clip located underneath the lid of the cartridge.
- Press down the cartridge to pierce the seal until fully engaged into the cartridge groove.
- Open the lid and check that the seal is completely punctured before loading a sample.
- Mix sample by rapidly inverting the specimen or external quality control tube 5 times. Carefully open the cap of the specimen tube or external quality control.
- Dispense 300 µL of the sample into the hole in the upper right corner of the cartridge using a 300 µL of fixed volume dropper or a micro pipette.


Figure 8. Fixed volume dropper

- After a few seconds, Sample Guide screen will automatically change to the Insert Cartridge screen. Touch the Sample Guide screen if you want to skip the guide.
- Close the lid.


Figure 9. Loading a sample

Figure 10. Sample Guide screen

Figure 11. Insert Cartridge screen

Running a test

1. Load the cartridge on the selected STANDARD M10 Module with the Amplification chamber facing the inside of the module. (The status indicator of the selected module will blink green.)
2. Close the door completely.
3. After confirming the sample and cartridge information, touch the OK button on the screen. (Touch the Reset button to re-input the information.)
4. Assay starts automatically, and remaining time will appear on the screen.

Figure 12. Confirm the test screen

Figure 13. Running screen

5. If the signal from the target reaches a certain threshold before the full process have been completed, an Early Detection Call function is initiated and it will provide earlier time to results (The Early Detection Call function is activated after 25 minutes of operation).

Figure 14. Early Detection Call function

6. When the run is finished, it switches to the Review screen and the result is displayed.
7. Dispose of used cartridges in the appropriate biohazard waste container according to your institution's standard practices.

8. To run another test, touch the Home icon and repeat the process (If another STANDARD M10 Module connected to the STANDARD M10 Console is available, you can start a new test while another test is running).

Interpretation of Results

The results are interpreted automatically by STANDARD M10 Console and are clearly shown in the Review screen. STANDARD M10 Flu/RSV/SARS-CoV-2 Fast test provides test results based on the detection of respective gene targets according to the algorithms shown in Table 3.

Table 3. Description of results

Outcome (Home screen)	Result (Review screen)	Description
Positive	+	At least one pathogen/target gene is positive.
Negative	—	No pathogen was detected.
Invalid	!	IC signal does not have a Ct within the valid range.
Error	X	The test failed because either an error occurred or the test was canceled by the user.

Table 4. Description of IC results

Outcome (Summary screen)	Result (Summary screen)	Description
IC Valid	V	IC has a Ct within the valid range. : The test was completed. Report positive/negative results of target according to the interpretation shown in table 5. (If the target is positive, it is displayed as IC Valid regardless of the IC Ct value.)
IC Invalid	!	IC does not meet acceptance criteria and all targets are not detected. Repeat test.
IC Error	X	The test failed because either an error occurred or the test was canceled by the user. Repeat the test.

Table 5. Interpretation of results

The STANDARD M10 Flu/RSV/SARS-CoV-2 Fast provides test results analyzed automatically by the STANDARD M10 Console using its proprietary algorithm.

Result	Interpretation
Flu A Positive	The Flu A target RNA is detected. - The Flu A signal has a Ct within the valid range. - IC: N/A (not applicable); IC may be inhibited and delayed Ct value when Flu A is positive at high concentration.
Flu B Positive	The Flu B target RNA is detected. - The Flu B signal has a Ct within the valid range. - IC: N/A (not applicable); IC may be inhibited and delayed Ct value when Flu B is positive at high concentration.
RSV Positive	The RSV target RNA is detected. - The RSV signal has a Ct within the valid range. - IC: N/A (not applicable); IC may be inhibited and delayed Ct value when RSV is positive at high concentration.
SARS-CoV-2 Positive	The SARS-CoV-2 target RNA is detected. - The SARS-CoV-2 signal has a Ct within the valid range. - IC: N/A (not applicable); IC may be inhibited and delayed Ct value when SARS-CoV-2 is positive at high concentration.
Flu A Positive, Flu B Positive, RSV Positive, SARS-CoV-2 Positive	- Flu A, Flu B, RSV or/and SARS-CoV-2 target RNAs are detected. - IC: N/A (not applicable); IC may be inhibited and delayed Ct value when Flu A, Flu B, RSV, and SARS-CoV-2 are positive at high concentration.
Flu A Negative, Flu B Negative, RSV Negative, SARS-CoV-2 Negative	- Flu A, Flu B, RSV and SARS-CoV-2 target RNAs are not detected. - IC: Valid; IC has a Ct within the valid range.
Invalid	IC does not meet acceptance criteria and all targets are not detected. Repeat test. IC: Invalid; IC and Flu A, Flu B, RSV and SARS-CoV-2 signals do not have a Ct within valid range.

Error	The test failed because either an error occurred or the test was canceled by the user. Presence or absence of target nucleic acids cannot be determined. Repeat the test.							
Interpretation	Influenza A/B			RSV		SARS-CoV-2		
	Flu A	Flu B	IC	RSV	IC	ORF1ab/E gene	N gene	IC
Flu A Positive	+	-	+/-	-	+	-	-	+
Flu B Positive	-	+	+/-	-	+	-	-	+
RSV Positive	-	-	+	+	+/-	-	-	+
SARS-CoV-2 Positive	-	-	+	-	+	+	+	+/-
	-	-	+	-	+	+	-	+/-
	-	-	+	-	+	-	+	+/-
Flu A Positive, Flu B Positive, RSV Positive, SARS-CoV-2 Positive	+	+	+/-	+	+/-	+	+	+/-
	+	+	+/-	+	+/-	+	-	+/-
	+	+	+/-	+	+/-	-	+	+/-
Flu A Negative, Flu B Negative, RSV Negative, SARS-CoV-2 Negative	-	-	+	-	+	-	-	+
Invalid	-	-	-	-	-	-	-	-

Figure 15. Individual result (left) and test history (right)

Note

- Co-infection can occur with these targets: Flu A, Flu B, RSV and SARS-CoV-2.
- This assay is not intended to differentiate Influenza A/B and RSV subtype or subgroups or lineage and SARS-CoV-2 variants. If differentiation of specific Influenza A/B or RSV or SARS-CoV-2 subtypes and strains is needed, additional testing, in consultation with state or local public health departments, is required.

Quality Control**1) Internal Control (IC)**

The internal control is intended to verify that the sample was properly applied, the cartridge reagents are functioning correctly, no interfering substances are present, and the procedure was performed accurately.

2) External Quality Control

External quality control procedure is intended to monitor cartridge and assay performance.

The specific test method for the external quality control is referred to the *Loading a sample into STANDARD M10 Flu/RSV/SARS-CoV-2 Fast cartridge* section. If the external quality control is not valid, the patient test result is not guaranteed and the result is invalid. If the expected results for the external control materials are not obtained, repeat the external controls, prior to releasing patient results. If the expected results for the external control material are not obtained upon repeat, contact Technical Support

The following external control is available but not provided. It may be used in accordance with the requirements of local, state, and federal accrediting organizations, as applicable:

- AccuPlex™ SARS-CoV-2, Flu A/B and RSV Molecular Controls Kit (SeraCare, 0505-0260)

Performance**■ Limit of Detection (LoD)**

The product was tested for Influenza A/B, RSV, and SARS-CoV-2 using 9 types of positive reference materials. These materials were serially diluted in negative nasopharyngeal swab samples and tested over 3 days using 2 different lots at 5 concentration levels, with 24 replicates per concentration. Based on the experimental results, the Limit of Detection (LoD) was determined through Probit analysis, and the LoD for each target is as follows.

Table 6. Limit of Detection (LoD) for STANDARD M10 Flu/RSV/SARS-CoV-2 Fast

Target Pathogen	Subtype (Strain)	Limit of Detection (95% CI)
Influenza A	H1N1pdm_A/NY/02/09	199.5 cp/mL (95% CI 117.5 - 331.1 cp/mL)
	H1_A/Singapore/63/04	354.8 cp/mL (95% CI 195.0 - 645.7 cp/mL)
	H3_A/Brisbane/10/07	141.3 cp/mL (95% CI 89.1 - 223.9 cp/mL)
Influenza B	Yamagata_Phuket/3073/13	0.013 TCID ₅₀ /mL (95% CI 0.008 - 0.019 TCID ₅₀ /mL)
	Victoria_Brisbane/33/08	0.016 TCID ₅₀ /mL (95% CI 0.009 - 0.028 TCID ₅₀ /mL)

RSV	RSV-A_SC2632	281.8 cp/mL (95% CI 154.9 - 501.2 cp/mL)	
	RSV-B_CH93(18)-18	104.7 cp/mL (95% CI 63.1 - 173.8 cp/mL)	
SARS-CoV-2	USA-WA1/2020	ORF1ab/E gene	263.0 cp/mL (95% CI 199.5 - 346.7 cp/mL)
		N gene	537.0 cp/mL (95% CI 380.2 - 758.6 cp/mL)
	WHO panel_England/02/2020	ORF1ab/E gene	182.0 IU/mL (95% CI 97.7 - 331.1 IU/mL)
		N gene	524.8 IU/mL (95% CI 281.8 - 1000.0 IU/mL)

■ Inclusivity

To evaluate the inclusivity of the product for its target viruses—Influenza A/B, RSV, and SARS-CoV-2—testing was conducted using 34 strains of Influenza A, 24 strains of Influenza B, 4 strains of RSV, and 31 variants of SARS-CoV-2. Each strain was diluted in negative nasopharyngeal swab samples to achieve a concentration of 1X LoD, based on the established Limit of Detection (LoD) for each target. The testing was performed with 5 replicates per strain, and the results confirmed the inclusivity of the product for the tested strains.

Table 7. Strains tested in Inclusivity test

Influenza A		Influenza B		RSV		SARS-CoV-2	
No.	Strain	No.	Strain	No.	Strain	No.	Strain
1	NY/03/09	1	Massachusetts/2/12	1	RSV A; Long	1	S. African B.1.351
2	New York/18/09	2	Wisconsin/1/10	2	RSV A; A2	2	Kappa B.1.617
3	A/California/07/2009	3	Florida/07/04	3	RSV A; 2006 Isolate	3	Delta B.1.617.2
4	Oseltamivir-Resistant (H275Y Isolate 2)	4	B/Maryland/1/59	4	RSV B; 18537	4	Delta Plus AY.1
5	Brisbane/2/18	5	B/Lee/40 TC adapted	-	-	5	Omicron B.1.1.529
6	Michigan/45/15	6	B/Florida/4/2006	-	-	6	Omicron BQ.1
7	Guangdong-Maonan/SWL 1536/19	7	Texas/6/11	-	-	7	Omicron BQ.1.1
8	Canada/6294/09	8	Yamagata/16/88	-	-	8	Omicron BA.4.6
9	New Caledonia/20/99	9	Utah/9/14	-	-	9	Omicron BA.5
10	Brisbane/59/07	10	Brisbane/9/14	-	-	10	Omicron BA.2.3
11	A/Virginia/ATCC1/2009	11	Brisbane/36/12	-	-	11	Omicron XBB

12	Solomon Islands/03/06	12	Panama/45/90	-	-	12	Omicron XBB.1.5
13	PR/8/34	13	B/GL/1739/54	-	-	13	Omicron BA.2.12.1
14	Taiwan/42/06	14	Malaysia/2506/04	-	-	14	Omicron BA.2.75.5
15	Texas/50/12	15	B/Florida/78/2015	-	-	15	Omicron BQ.1.3
16	Perth/16/09	16	Colorado/06/17	-	-	16	Omicron BA.5.5
17	A/Port Chalmers/1/73	17	Alabama/2/17	-	-	17	B.1
18	A/Victoria/3/75	18	Victoria/504/00	-	-	18	B.1.1.7
19	A/Aichi/2/68	19	Hawaii/10/01	-	-	19	B.1.526
20	HK/8/68	20	Brisbane/46/15	-	-	20	Zeta P2_2021
21	Wisconsin/67/05	21	Texas/2/13	-	-	21	Gamma P1_2021
22	California/7/04	22	Victoria/2/87	-	-	22	Omicron BF.5
23	Switzerland/9715293/13	23	Brisbane/33/08	-	-	23	Omicron BF.7
24	South Australia/55/14	24	B Austria /1359417/2021	-	-	24	Omicron EG.5.1
25	Stockholm/6/14	-	-	-	-	25	Omicron JN.1
26	Norway/466/14	-	-	-	-	26	Omicron JN.1.4
27	Kumamoto/102/02	-	-	-	-	27	Omicron JG.3
28	Singapore/ INFIMH-16-0019/16	-	-	-	-	28	Omicron CA.1
29	Kansas/14/17	-	-	-	-	29	Omicron BA.5.2.6
30	Victoria/361/11	-	-	-	-	30	Omicron BQ.1.16
31	Texas/50/12	-	-	-	-	31	B.1.595
32	Brisbane/10/07	-	-	-	-	-	-
33	Hong Kong/4801/14	-	-	-	-	-	-
34	A/Darwin/6/2021	-	-	-	-	-	-

■ Cross-reactivity & Microbial Interference

A total of 46 cross-reactive substances, including those with genetic relevance to the product's targets or those causing similar symptoms, were selected for testing. These substances were diluted in negative nasopharyngeal swab samples to the test concentration and evaluated in three replicates at both negative sample levels and 3X LoD concentrations for each target. The results confirmed that there were no cross-reactivity or microbial interferences at the tested concentrations for the 46 substances listed below.

Table 8. Substances tested in Cross-reactivity & Microbial Interference test

No.	Cross-reactive (microbial interference) Substances	Test Concentration	No.	Cross-reactive (microbial interference) Substances	Test Concentration
1	MERS-CoV	1.00x10 ⁵ TCID ₅₀ /mL	24	<i>Moraxella catarrhalis</i>	5.92x10 ⁶ CFU/mL
2	HCoV-OC43	1.00x10 ⁵ TCID ₅₀ /mL	25	Cytomegalovirus (CMV)	1.00x10 ⁶ cp/mL
3	HCoV-NL63	1.00x10 ⁵ TCID ₅₀ /mL	26	Measles virus	1.00x10 ⁵ TCID ₅₀ /mL
4	HCoV-229E	1.00x10 ⁵ TCID ₅₀ /mL	27	Mumps virus	1.15x10 ⁶ TCID ₅₀ /mL
5	Adenovirus Type03 (Species B)	1.00x10 ⁵ TCID ₅₀ /mL	28	Epstein Barr virus	1.00x10 ⁶ cp/mL
6	Adenovirus Type01 (Species C)	3.72x10 ⁶ TCID ₅₀ /mL	29	<i>Candida albicans</i>	1.00x10 ⁶ CFU/mL
7	Human Metapneumovirus (hMPV)	1.00x10 ⁵ TCID ₅₀ /mL	30	<i>Legionella pneumophila</i>	1.00x10 ⁷ CFU/mL
8	Parainfluenza Virus Type 1	1.00x10 ⁵ TCID ₅₀ /mL	31	<i>Mycobacterium tuberculosis</i>	1.00x10 ⁶ cp/mL
9	Parainfluenza Virus Type 2	1.00x10 ⁵ TCID ₅₀ /mL	32	<i>Bordetella parapertussis</i>	1.00x10 ⁷ CFU/mL
10	Parainfluenza Virus Type 3	3.39x10 ⁵ TCID ₅₀ /mL	33	<i>Aspergillus sp.</i>	1.00x10 ⁶ CFU/mL
11	Parainfluenza Virus Type 4A	3.39x10 ⁵ TCID ₅₀ /mL	34	<i>Fusobacterium necrophorum</i>	1.00x10 ⁶ CFU/mL
12	Parainfluenza Virus Type 4B	1.00x10 ⁵ TCID ₅₀ /mL	35	<i>Corynebacterium sp.</i>	3.00x10 ⁶ CFU/mL
13	Enterovirus Type 68	1.00x10 ⁵ TCID ₅₀ /mL	36	<i>Escherichia coli</i>	1.00x10 ⁷ CFU/mL
14	Rhinovirus A16	1.00x10 ⁵ TCID ₅₀ /mL	37	<i>Lactobacillus sp.</i>	1.00x10 ⁷ CFU/mL
15	Rhinovirus B	4.00x10 ⁵ TCID ₅₀ /mL	38	<i>Neisseria meningitidis</i>	1.00x10 ⁶ CFU/mL
16	Haemophilus Influenzae Type B	1.42x10 ⁶ CFU/mL	39	<i>Neisseria gonorrhoeae</i> , titrated (Strain : 2017)	1.00x10 ⁶ CFU/mL
17	<i>Streptococcus pneumoniae</i>	1.42x10 ⁶ CFU/mL	40	<i>Staphylococcus aureus</i>	1.00x10 ⁷ CFU/mL
18	<i>Streptococcus pyogenes</i>	4.60x10 ⁶ CFU/mL	41	<i>Mycoplasma genitalium</i>	1.00x10 ⁶ CFU/mL
19	<i>Bordetella pertussis</i>	2.21x10 ⁷ CFU/mL	42	<i>Chlamydomphila pneumoniae</i>	1.00x10 ⁶ IFU/mL
20	<i>Mycoplasma pneumoniae</i>	1.40x10 ⁶ CFU/mL	43	<i>Pneumocystis jirovecii</i> (PJP)	1.00x10 ⁷ CFU/mL

21	<i>Pseudomonas aeruginosa</i>	5.60x10 ⁶ CFU/mL	44	Human coronavirus HKU1	1.00x10 ⁶ cp/mL
22	<i>Staphylococcus epidermis</i>	1.20x10 ⁶ CFU/mL	45	Pooled human nasal wash	Stock
23	<i>Streptococcus salivarius</i>	1.35x10 ⁶ CFU/mL	46	SARS-CoV	*1.00x10 ⁵ cp/mL
-	-	-			**1.00x10 ⁶ cp/mL

* Cross-reactivity with SARS-CoV-2 N gene, but no microbial interference reaction

** No cross-reactivity with SARS-CoV-2 N gene

■ Interfering Substance

A total of 17 endogenous and exogenous interfering substances, which may potentially be present in the specimen, were tested. These substances were diluted to the test concentration using negative nasopharyngeal swab samples and evaluated under conditions with and without the addition of interfering substances in both negative and 3X LoD concentration samples, with three replicates per condition. The results confirmed that there was no interference up to the specified test concentrations for the 17 substances listed below.

Table 9. Interfering substances tested in Interfering Substance test

No.	Interfering Substance	Active Ingredient	Test Concentration
1	Human whole Blood	Plasma, Blood cell	5% v/v
2	Human genomic DNA	human genome, protein	15 µg/mL
3	Hemoglobin	Hemoglobin	0.2 g/mL
4	Mucin	Mucin	2.5 mg/mL
5	Albumin	Albumin	0.24 mg/mL
6	Biotin	Biotin	200 µg/mL
7	Dexamethasone	Dexamethasone	6 mg/mL
8	Mouthwash	Cetylpyridinium chloride	5% v/v
9	Oseltamivir	Oseltamivir	25 mg/mL
10	Oxymetazoline (Nasal Spray)	Oxymetazoline	15% v/v
11	Ribavirin	Ribavirin	800 mg/mL
12	Chloraseptic Spray (Phenol/Glycerin)	Phenol	15% v/v
13	Sore throat and cough lozenges	Benzocaine	3 mg/mL
14	Tobacco	Nicotine, Tar	0.03 mg/mL

15	Tobramycin	Tobramycin	512 µg/mL
16	Toothpaste	Sodium monofluorophosphate	0.5% v/v
17	Zanamivir	Zanamivir	3.3 mg/mL

■ Competitive Interference

To evaluate the competitive interference response between the analytes, the positive reference materials for Influenza A/B, RSV, and SARS-CoV-2 targets of low concentration (3X LoD) and high concentration (100X LoD) were prepared and evaluated. Based on the test results, there was no competitive interference reaction between the analytes.

■ Precision (Repeatability & Reproducibility)

[Repeatability]

The positive reference materials for Influenza A/B, RSV, and SARS-CoV-2 targets were diluted in negative nasopharyngeal swab samples to three test concentrations (0, 1X, 3X LoD). Testing was performed with one lot, four repetitions per day, over twenty days. Based on the test results, SD values for within-run, between-run, between-day, and within-laboratory were analyzed, confirming the repeatability of the product with an SD < 2.0 Ct.

[Reproducibility]

The positive reference materials for Influenza A/B, RSV, and SARS-CoV-2 targets were diluted in negative nasopharyngeal swab samples to three test concentrations (0, 1X, 3X LoD). Testing was performed with between-operator, between-lots, between-site, and instrument five repetitions per day over five days. Based on the test results, the SD and CV (%) values were analyzed, confirming the reproducibility of the product with an SD < 2.0 Ct and CV < 5%.

■ Carry-over / Cross-contamination

The high-concentration positive reference materials for Influenza A/B, RSV, and SARS-CoV-2 were tested in cross-validation using negative nasopharyngeal swab samples with one lot of the product, repeating the test five times. The results showed that all high-concentration positive samples were positive, and all negative samples were negative, confirming that there was no cross-contamination between the dedicated equipment used to run the product (STANDARD M10 System) and the test product.

■ Clinical Performance Study

The clinical performance of the M10 Flu/RSV/SARS-CoV-2 Fast test was evaluated using archive clinical nasopharyngeal (NP) swab specimens in Copan UTM. A total of 722 NP swab specimens were tested with M10 Flu/RSV/SARS-CoV-2 Fast by comparing its results with those from a state-of-the-art, CE-marked RT-PCR comparator device.

M10 Flu/RSV/SARS-CoV-2 Fast demonstrated a sensitivity and specificity of 100.00% and 99.13% for Influenza A, respectively; 98.78% and 98.44% for Influenza B, respectively; 97.96% and 98.72% for RSV, respectively; 98.51% and 99.49% for SARS-CoV-2, respectively (Table 10).

Table 10. Summary of the clinical sensitivity and specificity test results

1) Influenza A virus

- Clinical sensitivity: 100.00% (148/148) [95% CI: 97.47% - 100.00%]

- Clinical specificity: 99.13% (569/574) [95% CI: 97.98% - 99.63%]

		Confirmed result		Total
		Positive	Negative	
STANDARD M10 Flu/RSV/SARS-CoV-2 Fast	Positive	148	5	153
	Negative	0	569	569
	Total	148	574	722

2) Influenza B virus

- Clinical sensitivity: 98.78% (81/82) [95% CI: 93.41% - 99.78%]

- Clinical specificity: 98.44% (630/640) [95% CI: 97.15% - 99.15%]

		Confirmed result		Total
		Positive	Negative	
STANDARD M10 Flu/RSV/SARS-CoV-2 Fast	Positive	81	10	91
	Negative	1	630	631
	Total	82	640	722

3) Respiratory Syncytial virus

- Clinical sensitivity: 97.96% (96/98) [95% CI: 92.86% - 99.44%]

- Clinical specificity: 98.72% (616/624) [95% CI: 97.49% - 99.35%]

		Confirmed result		Total
		Positive	Negative	
STANDARD M10 Flu/RSV/SARS-CoV-2 Fast	Positive	96	8	104
	Negative	2	616	618
	Total	98	624	722

4) SARS-CoV-2

- Clinical sensitivity: 98.51% (132/134) [95% CI: 94.72% - 99.59%]

- Clinical specificity: 99.49% (585/588) [95% CI: 98.51% - 99.83%]

		Confirmed result		Total
		Positive	Negative	
STANDARD M10 Flu/RSV/SARS-CoV-2 Fast	Positive	132	3	135
	Negative	2	585	587
	Total	134	588	722

[Clinical Performance in Near-Patient Testing (NPT) Setting]

The clinical performance of the STANDARD™ M10 Flu/RSV/SARS-CoV-2 Fast was evaluated in a near-patient testing (NPT) environment using nasopharyngeal (NP) swab specimens.

The performance in the NPT setting was compared to the performance in the laboratory environment by clinical laboratory professionals or healthcare professionals who have received general professional training appropriate to their role and the usability in the NPT setting by intended users (HCPs).

A total of 379 NP swab samples were analyzed. Each sample was tested with the M10 test in both NPT and laboratory settings, and results were compared with a CE-marked RT-PCR reference test and resolver assay.

Table 11. Summary of the clinical performance results in NPT setting

1) Influenza A Positive (Total sample : 96)

- NPT setting compared to private lab setting

Influenza A target		Private Lab Setting	
		Influenza A Positive	Negative
NPT-setting	Influenza A Positive	96	0
	Negative	0	0

Positive percentage agreement (PPA) : 100% (Wilson 95% CI : 95.36-100%)

- NPT setting compared to reference/resolver assay

Influenza A target		Reference/resolver Assay	
		Influenza A Positive	Negative
NPT-setting	Influenza A Positive	96	0
	Negative	0	0

Diagnostic sensitivity : 100% (Wilson 95% CI : 95.36-100%)

- 2) Influenza B Positive (Total sample : 54)
- NPT setting compared to private lab setting

Influenza B target		Private Lab Setting	
		Influenza B Positive	Negative
NPT-setting	Influenza B Positive	52	0
	Negative	2	0
Positive percentage agreement (PPA) : 96.15% (Wilson 95% CI : 87.02-98.94%)			

- NPT setting compared to reference/resolver assay

Influenza B target		Reference/resolver Assay	
		Influenza B Positive	Negative
NPT-setting	Influenza B Positive	52	0
	Negative	2	0
Diagnostic sensitivity : 96.15% (Wilson 95% CI : 87.02-98.94%)			

- 3) RSV target Positive (Total sample : 64)
- NPT setting compared to private lab setting

RSV target		Private Lab Setting	
		RSV Positive	Negative
NPT-setting	RSV Positive	61	3
	Negative	0	0
Positive percentage agreement (PPA) : 95.31% (Wilson 95% CI : 87.10-98.39%)			

- NPT setting compared to reference/resolver assay

RSV target		Reference/resolver Assay	
		RSV Positive	Negative
NPT-setting	RSV Positive	64	0
	Negative	0	0
Diagnostic sensitivity : 100% (Wilson 95% CI : 95.34-100%)			

- 4) SARS-CoV-2 Positive (Total sample : 63)

- NPT setting compared to private lab setting

SARS-CoV-2 target		Private Lab Setting	
		SARS-CoV-2 Positive	Negative
NPT-setting	SARS-CoV-2 Positive	63	0
	Negative	0	0
Positive percentage agreement (PPA) : 100% (Wilson 95% CI : 94.25-100%)			

- NPT setting compared to reference/resolver assay

SARS-CoV-2 target		Reference/resolver Assay	
		SARS-CoV-2 Positive	Negative
NPT-setting	SARS-CoV-2 Positive	63	0
	Negative	0	0
Positive percentage agreement (PPA) : 100% (Wilson 95% CI : 94.25-100%)			

- 5) Negative Sample (Total sample : 119)

- NPT setting compared to private lab setting

All the targets		Private Lab Setting	
		Positive	Negative
NPT-setting	Positive	0	2
	Negative	2	115
Negative percentage agreement (NPA) : 96.64% (Wilson 95% CI : 91.68-98.69%)			

- NPT setting compared to reference/resolver assay

All the targets		Reference/Resolver Assay	
		Positive	Negative
NPT-setting	Positive	0	2
	Negative	2	117
Diagnostic specificity : 98.32% (Wilson 95% CI : 94.08-99.54%)			

Limitations

1. This product must be used with human nasopharyngeal swab specimens.
2. False negative results may occur from several causes.
 - Sample concentrations is near or below the limit of detection of the test.
 - A specimen is improperly collected, transported or handled.
 - Cartridges are exposed to improper environmental factors.
3. Modifications to the procedures provided in the IFU (Instructions for Use) within the package may alter the performance of the test.
4. False positive results may occur from cross-contamination between patient samples, specimen mix-up and/or RNA contamination during product handling.
5. Qualitative detection of positive results in this cartridge does not indicate the presence of infectious virus.
6. This cartridge only classifies and identifies the Flu A, Flu B, RSV or SARS-CoV-2. The test results are for clinical reference only, and it should not be used as the only evidence for clinical diagnosis and treatment.
7. Mutations within the target regions of the M10 Flu/RSV/SARS-CoV-2 Fast test could affect primer and/or probe binding resulting in failure to detect the presence of virus or the virus being detected less predictably.
8. This product has the potential to cross-react with SARS-CoV-1 in the N gene of SARS-CoV-2. However, no cases of infection with this pathogen have been reported since 2004 (Source: CDC), and human infections are extremely rare. Therefore, the risk associated with the test results is considered low.

References

1. Characteristics of SARS-CoV-2 and COVID-19. *Nature Reviews Micro.*, 2020, <http://doi.org/10.1038/s41579-020-00459-7>
2. Diagnosis and treatment of pneumonia caused by new coronavirus (trial version 4) National Health Commission. 2020
3. Centers for Disease Control and Prevention. 2020. CDC guidance for expanded screening testing to reduce silent spread of SARS-CoV-2. CDC, Atlanta, GA. https://www.cdc.gov/csels/dls/locs/2020/cdc_publishes_new_guidance_for_expanded_screening_testing_to_reduce_silent_spread_of_sars-cov-2.html.
4. Centers for Disease Control and Prevention. 2020. Overview of testing for SARS-CoV-2 (COVID-19). <https://www.cdc.gov/coronavirus/2019-ncov/hcp/testing-overview.html>. Accessed 3 June 2021.
5. Johns Hopkins University of Medicine. 2021. Coronavirus Resource Center. <https://coronavirus.jhu.edu/map.html>. Accessed 3 June 2021.
6. Belongia EA, Osterholm MT. 2020. COVID-19 and flu, a perfect storm. *Science* 368:1163. <https://doi.org/10.1126/science.abd2220>.
7. Zhang N, Wang L, Deng X, Liang R, Su M, He C, Hu L, Su Y, Ren J, Yu F, Du L, Jiang S. 2020. Recent advances in the detection of respiratory virus infection in humans. *J Med Virol* 92:408–417. <https://doi.org/10.1002/jmv.25674>.
8. Jiang C, Yao X, Zhao Y, Wu J, Huang P, Pan C, Liu S, Pan C. 2020. Comparative review of respiratory diseases caused by coronaviruses and influenza A viruses during epidemic season. *Microbes Infect* 22:236–244. <https://doi.org/10.1016/j.micinf.2020.05.005>.

9. Obermeier M, Pacenti M, Ehret R, Onelia F, Gunson R, Goldstein E, Chevaliez S, Vilas A, Glass A, Maree L, Krügel M, Knechten H, Braun P, Naeth G, Azzato F, Lucic D, Marlowe N, Palm MJ, Pfeifer K, Reinhardt B, Dhein J, Joseph AM, Martínez-García L, Galán J-C. 2020. Improved molecular laboratory productivity by consolidation of testing on the new random-access analyzer Alinity m. *J Lab Med* 44:319–328. <https://doi.org/10.1515/labmed-2020-0102>.
10. Walsh EE. 2017. Respiratory syncytial virus infection: an illness for all ages. *Clin Chest Med* 38:29–36. <https://doi.org/10.1016/j.ccm.2016.11.010>.
11. REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on the classification labeling and packaging of substances and mixtures amending and repealing, List of Precautionary Statements, Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

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Omezení

1. This product must be used with human nasopharyngeal swab specimens.
2. False Negativní results may occur from several causes.
 - Sample concentrations is near or below the limit of detection of the test.
 - A specimen is improperly collected, transported or handled.
 - Cartridges are exposed to improper environmental factors.
3. Modifications to the procedures provided in the IFU (Instructions for Use) within the package may alter the performance of the test.
4. False Pozitivní results may occur from cross-contamination between patient samples, specimen mix-up and/or RNA contamination during product handling.
5. Qualitative detection of Pozitivní results in this cartridge does not indicate the presence of infectious virus.
6. This cartridge only classifies and identifies the Flu A, Flu B, RSV or SARS-CoV-2. The test results are for clinical reference only, and it should not be used as the only evidence for clinical diagnosis and treatment.
7. Mutations within the target regions of the M10 Flu/RSV/SARS-CoV-2 Fast test could affect primer and/or probe binding resulting in failure to detect the presence of virus or the virus being detected less predictably.
8. This product has the potential to cross-react with SARS-CoV-1 in the N gene of SARS-CoV-2. However, no cases of infection with this pathogen have been reported since 2004 (Source: CDC), and human infections are extremely rare. Therefore, the risk associated with the test results is considered low.

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1. Characteristics of SARS-CoV-2 and COVID-19. Nature Reviews Micro., 2020, <http://doi.org/10.1038/s41579-020-00459-7>
2. Diagnosis and treatment of pneumonia caused by new coronavirus (trial version 4) National Health Commission. 2020
3. Centers for Disease Control and Prevention. 2020. CDC guidance for expanded screening testing to reduce silent spread of SARS-CoV-2. CDC, Atlanta, GA. https://www.cdc.gov/csels/dls/locs/2020/cdc_publicshes_new_guidance_for_expanded_screening_testing_to_reduce_silent_spread_of_sars-cov-2.html.
4. Centers for Disease Control and Prevention. 2020. Overview of testing for SARS-CoV-2 (COVID-19). <https://www.cdc.gov/coronavirus/2019-ncov/hcp/testing-overview.html>. Accessed 3 June 2021.
5. Johns Hopkins University of Medicine. 2021. Coronavirus Resource Center. <https://coronavirus.jhu.edu/map.html>. Accessed 3 June 2021.
6. Belongia EA, Osterholm MT. 2020. COVID-19 and flu, a perfect storm. Science 368:1163. <https://doi.org/10.1126/science.abd2220>.
7. Zhang N, Wang L, Deng X, Liang R, Su M, He C, Hu L, Su Y, Ren J, Yu F, Du L, Jiang S. 2020. Recent advances in the detection of respiratory virus infection in humans. J Med Virol 92:408–417. <https://doi.org/10.1002/jmv.25674>.
8. Jiang C, Yao X, Zhao Y, Wu J, Huang P, Pan C, Liu S, Pan C. 2020. Comparative review of respiratory diseases caused by coronaviruses and influenza A viruses during epidemic season. Microbes Infect 22:236–244. <https://doi.org/10.1016/j.micinf.2020.05.005>.
9. Obermeier M, Pacenti M, Ehret R, Onelia F, Gunson R, Goldstein E, Chevaliez S, Vilas A, Glass A, Maree L, Krügel M, Knecht H, Braun P, Naeth G, Azzato F, Lucic D, Marlowe N, Palm MJ, Pfeifer K, Reinhardt B, Dhein J, Joseph AM, Martínez-García L, Galán J-C. 2020. Improved molecular laboratory productivity by consolidation of testing on the new random-access analyzer Alinity m. J Lab Med 44:319–328. <https://doi.org/10.1515/labmed-2020-0102>.
10. Walsh EE. 2017. Respiratory syncytial virus infection: an illness for all ages. Clin Chest Med 38:29–36. <https://doi.org/10.1016/j.ccm.2016.11.010>.
11. REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on the classification labeling and packaging of substances and mixtures amending and repealing, List of Precautionary Statements, Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

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Referencie

- Characteristics of SARS-CoV-2 and COVID-19. Nature Reviews Micro., 2020, <http://doi.org/10.1038/s41579-020-00459-7>
- Diagnosis and treatment of pneumonia caused by new coronavirus (trial version 4) National Health Commission. 2020

- Centers for Disease Control and Prevention. 2020. CDC guidance for expanded screening testing to reduce silent spread of SARS-CoV-2. CDC, Atlanta, GA. https://www.cdc.gov/csels/dls/locs/2020/cdc_publishes_new_guidance_for_expanded_screening_testing_to_reduce_silent_spread_of_sars-cov-2.html.
- Centers for Disease Control and Prevention. 2020. Overview of testing for SARS-CoV-2 (COVID-19). <https://www.cdc.gov/coronavirus/2019-ncov/hcp/testing-overview.html>. Accessed 3 June 2021.
- Johns Hopkins University of Medicine. 2021. Coronavirus Resource Center. <https://coronavirus.jhu.edu/map.html>. Accessed 3 June 2021.
- Belongia EA, Osterholm MT. 2020. COVID-19 and flu, a perfect storm. Science 368:1163. <https://doi.org/10.1126/science.abd2220>.
- Zhang N, Wang L, Deng X, Liang R, Su M, He C, Hu L, Su Y, Ren J, Yu F, Du L, Jiang S. 2020. Recent advances in the detection of respiratory virus infection in humans. J Med Virol 92:408–417. <https://doi.org/10.1002/jmv.25674>.
- Jiang C, Yao X, Zhao Y, Wu J, Huang P, Pan C, Liu S, Pan C. 2020. Comparative review of respiratory diseases caused by coronaviruses and influenza A viruses during epidemic season. Microbes Infect 22:236–244. <https://doi.org/10.1016/j.micinf.2020.05.005>.
- Obermeier M, Pacenti M, Ehret R, Onelia F, Gunson R, Goldstein E, Chevaliez S, Vilas A, Glass A, Maree L, Krügel M, Knechten H, Braun P, Naeth G, Azzato F, Lucic D, Marlowe N, Palm MJ, Pfeifer K, Reinhardt B, Dhein J, Joseph AM, Martínez-García L, Galán J-C. 2020. Improved molecular laboratory productivity by consolidation of testing on the new random-access analyzer Alinity m. J Lab Med 44:319–328. <https://doi.org/10.1515/labmed-2020-0102>.
- Walsh EE. 2017. Respiratory syncytial virus infection: an illness for all ages. Clin Chest Med 38:29–36. <https://doi.org/10.1016/j.ccm.2016.11.010>.
- REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on the classification labeling and packaging of substances and mixtures amending and repealing, List of Precautionary Statements, Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

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Bibliografia

1. Characteristics of SARS-CoV-2 and COVID-19. *Nature Reviews Micro.*, 2020, <http://doi.org/10.1038/s41579-020-00459-7>
2. Diagnosis and treatment of pneumonia caused by new coronavirus (trial version 4) National Health Commission. 2020
3. Centers for Disease Control and Prevention. 2020. CDC guidance for expanded screening testing to reduce silent spread of SARS-CoV-2. CDC, Atlanta, GA. https://www.cdc.gov/csels/dls/locs/2020/cdc_publishes_new_guidance_for_expanded_screening_testing_to_reduce_silent_spread_of_sars-cov-2.html.
4. Centers for Disease Control and Prevention. 2020. Overview of testing for SARS-CoV-2 (COVID-19). <https://www.cdc.gov/coronavirus/2019-ncov/hcp/testing-overview.html>. Accessed 3 June 2021.
5. Johns Hopkins University of Medicine. 2021. Coronavirus Resource Center. <https://coronavirus.jhu.edu/map.html>. Accessed 3 June 2021.
6. Belongia EA, Osterholm MT. 2020. COVID-19 and flu, a perfect storm. *Science* 368:1163. <https://doi.org/10.1126/science.abd2220>.
7. Zhang N, Wang L, Deng X, Liang R, Su M, He C, Hu L, Su Y, Ren J, Yu F, Du L, Jiang S. 2020. Recent advances in the detection of respiratory virus infection in humans. *J Med Virol* 92:408–417. <https://doi.org/10.1002/jmv.25674>.
8. Jiang C, Yao X, Zhao Y, Wu J, Huang P, Pan C, Liu S, Pan C. 2020. Comparative review of respiratory diseases

caused by coronaviruses and influenza A viruses during epidemic season. *Microbes Infect* 22:236–244. <https://doi.org/10.1016/j.micinf.2020.05.005>.

9. Obermeier M, Pacenti M, Ehret R, Onelia F, Gunson R, Goldstein E, Chevaliez S, Vilas A, Glass A, Maree L, Krügel M, Knechten H, Braun P, Naeth G, Azzato F, Lucic D, Marlowe N, Palm MJ, Pfeifer K, Reinhardt B, Dhein J, Joseph AM, Martínez-García L, Galán J-C. 2020. Improved molecular laboratory productivity by consolidation of testing on the new random-access analyzer Alinity m. *J Lab Med* 44:319–328. <https://doi.org/10.1515/labmed-2020-0102>.
10. Walsh EE. 2017. Respiratory syncytial virus infection: an illness for all ages. *Clin Chest Med* 38:29–36. <https://doi.org/10.1016/j.ccm.2016.11.010>.
11. REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on the classification labeling and packaging of substances and mixtures amending and repealing, List of Precautionary Statements, Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

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Referanslar


1. Characteristics of SARS-CoV-2 and COVID-19. Nature Reviews Micro., 2020, <http://doi.org/10.1038/s41579-020-00459-7>
2. Diagnosis and treatment of pneumonia caused by new coronavirus (trial version 4) National Health Commission. 2020
3. Centers for Disease Control and Prevention. 2020. CDC guidance for expanded screening testing to reduce silent spread of SARS-CoV-2. CDC, Atlanta, GA. https://www.cdc.gov/csels/dls/locs/2020/cdc_publishes_new_guidance_for_expanded_screening_testing_to_reduce_silent_spread_of_sars-cov-2.html.
4. Centers for Disease Control and Prevention. 2020. Overview of testing for SARS-CoV-2 (COVID-19). <https://www.cdc.gov/coronavirus/2019-ncov/hcp/testing-overview.html>. Accessed 3 June 2021.
5. Johns Hopkins University of Medicine. 2021. Coronavirus Resource Center. <https://coronavirus.jhu.edu/map.html>. Accessed 3 June 2021.
6. Belongia EA, Osterholm MT. 2020. COVID-19 and flu, a perfect storm. Science 368:1163. <https://doi.org/10.1126/science.abd2220>.
7. Zhang N, Wang L, Deng X, Liang R, Su M, He C, Hu L, Su Y, Ren J, Yu F, Du L, Jiang S. 2020. Recent advances in the detection of respiratory virus infection in humans. J Med Virol 92:408–417. <https://doi.org/10.1002/jmv.25674>.
8. Jiang C, Yao X, Zhao Y, Wu J, Huang P, Pan C, Liu S, Pan C. 2020. Comparative review of respiratory diseases caused by coronaviruses and influenza A viruses during epidemic season. Microbes Infect 22:236–244. <https://doi.org/10.1016/j.micinf.2020.05.005>.
9. Obermeier M, Pacenti M, Ehret R, Onelia F, Gunson R, Goldstein E, Chevaliez S, Vilas A, Glass A, Maree L, Krügel M, Knechten H, Braun P, Naeth G, Azzato F, Lucic D, Marlowe N, Palm MJ, Pfeifer K, Reinhardt B, Dhein J, Joseph AM, Martinez-García L, Galán J-C. 2020. Improved molecular laboratory productivity by consolidation of testing on the new random-access analyzer Alinity m. J Lab Med 44:319–328. <https://doi.org/10.1515/labmed-2020-0102>.
10. Walsh EE. 2017. Respiratory syncytial virus infection: an illness for all ages. Clin Chest Med 38:29–36. <https://doi.org/10.1016/j.ccm.2016.11.010>.
11. REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on the classification labeling and packaging of substances and mixtures amending and repealing, List of Precautionary Statements, Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

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SD Biosensor, Inc.

C-4th&5th, 16, Deogyong-daero 1556beon-gil, Yeongtong-gu, Suwon-si, Gyeonggi-do, 16690,
REPUBLIC OF KOREA

Manufacturing site

14, Jeungpyeongsandan-ro, Jeungpyeong-eup, Jeungpyeong-gun, Chungcheongbuk-do, 27915,
REPUBLIC OF KOREA

ts@sdbiosensor.com | +82-31-300-0400 | www.sdbiosensor.com

MT Promed+Consulting GmbH

Ernst-Heckel-Straße 7 66386 St. Ingbert Germany
+49 6894 581020 | Fax: +49 6894 581021