# BAMBANKER® Cell Freezing Medium (a patented product)



#### Advantages of serum-free media

- · Stable quality. Composition-related variations between lots are limited as compared to serum-containing media.
- No risk of exposure of cells to unknown biological components or infectious substances.
- Serum-starved cells can be frozen. No reconditioning is needed.


Storage: Cool, dark place (at 2–10°C)

Shelf life: 24 months from the date of manufacture

#### Precautions

- Use cells during the logarithmic growth phase.
- Cannot be used for the human body.
- Do not use for any purpose other than research.
- Contains compounds derived from bovine serum.
- Before using this product, conduct a confirmation test using the desired cells.
- · We shall not be responsible for accidents or damage resulting from use of this product.

# **Cell freezing method using BAMBANKER**



Transfer the cell suspension to a centrifuging tube  $(5 \times 10^{5} 1 \times 10^7$  cells)

Remove the supernatant

Suspend the cells in 1 mL of suspension BAMBANKER into cryotube

## **Applications**

•3T3 - L1 •KLM - 1

-A431 ·MDCK

BAEC MEF

·Balb/3T3 ·NIH3T3

-C2C12 ·OKT4

•OP9 Daudi

•ECV304 •P3U1

PANC - 1 •H295R

·HEK293 •PC12

·HEK293T ·RPE

·HeLa ·SNL

·HeLa S3 ·TSU - Pr1

·HepG2 Vero

·HFF ·Human γ δ T cells

·Huh7 - Activated human peripheral blood lymphocytes

 Jurkat - Immortalized human muscle cells

·K562 · Mouse ES cell lines

· KATOIII Activated mouse spleen lymphocytes

#### This product has been successfully used with ES cells (mouse).

T.Hikichi, et al; Differentiation Potential of Parthenogenetic Embryonic

Stem Cells Is Improved

by Nuclear Transfer, Stem Cells, 2007, 25, 46-53

Code No.	Manufacturer code	Name	Quantity
302-14681	CS-02-001	BAMBANKER®	120 mL
302-14684	CS-04-001		20 mL x 5 bottles

Aliquot the

# **BAMBANKER®** Direct Cell Freezing Medium



## Advantages of serum-free media

- Stable quality. Composition-related variations between lots are limited as compared to serum-containing media.
- No risk of exposure of cells to unknown biological components or infectious substances.
- Serum-starved cells can be frozen. No reconditioning is needed.

Storage: Cool, dark place (at 2–10°C)


Shelf life: 24 months from the date of manufacture

#### Precautions


- Use cells during the logarithmic growth phase.
- Cannot be used for the human body.
- Do not use for any purpose other than research.
- Contains compounds derived from bovine serum.
- · Before using this product, conduct a confirmation test using the desired cells
- We shall not be responsible for accidents or damage resulting from use of this product.

## Cell freezing method using BAMBANKER® Direct

Ordinary preservation media







amount of
BAMBANKER
Direct as the
culture medium

Aliquot the suspension into cryotube

Code No.	Manufacturer code	Name	Quantity
306-95921	CS-06-001	BAMBANKER® Direct	20 mL

# Wako

#### Wako Pure Chemical Industries, Ltd.

Head Office: 1-2 Doshomachi 3-chome, Chuo-ku, Osaka 540-8605, Japan

Tokyo Main Office: 4-1 Nihonbashi Honcho 2-chome, Chuo-ku, Tokyo 103-0023, Japan

Phone 0120-052-099, fax 0120-052-806 (toll-free)

URL: http://www.wako-chem.co.jp

## Manufacturer



Lymphotec Inc.

18-4 Fuyuki, Koto-ku, Tokyo 135-0041, Japan

Phone: 03-3630-2530 (main) URL: http://www.lymphotec.co.jp