3diag - sTfR - CONTROL

REF TD-42683

LOT **683-08**

🕒 2023-03

(01) 08434477510052
 (10) 683-08
 (17) 230330
 (240) TD-42683

Calibrator & Method	Level	REF	mg/l target (range)	nmol/l target (range)
using 3diag - sTfR - CAL SET NIPIA, on IMMAGE®, and PETIA (R2 as starter)	CONTROL H sTfR ●	TD-42683-H	3.09 (2.47 - 3.71)	43.2 (34.5 - 51.8)
	CONTROL L sTfR ○	TD-42683-L	1.63 (1.30 - 1.95)	22.7 (18.2 - 27.2)
using 3diag - sTfR - CAL PETIA (sample as starter)	CONTROL H sTfR ●	TD-42683-H	2.96 (2.37 - 3.55)	41.3 (33.1 - 49.6)
	CONTROL L sTfR ○	TD-42683-L	1.37 (1.10 - 1.64)	19.1 (15.3 - 22.9)

- (EN)** Soluble Transferrin Receptor - Controls (2 lev.)
- (ES)** Receptor Soluble de la Transferrina - Controles (2 niv.)
- (IT)** Recettore Solubile della Transferrina - Controlli (2 liv.)
- (NL)** Oplosbare Transferrinereceptor - Controles (2 niv.)
- (PT)** Receptor Solúvel da Transferrina - Controlos (2 niv.)

TABLE of ANALYTICAL VALUES
TABLA de VALORES ANALITICOS
TABELLA di VALORI ANALITICI
ANALYTISCHE WAARDENTABEL
TABELA de VALORES ANALÍTICOS

<https://www.3diag.com/510>
 for IFU, Scan or follow link,
 and select LOT

IMMAGE® is a registered trademark of Beckman Coulter, Inc., Fullerton, CA

CONTROLS INCLUDED IN  **Accuracy 365** Online Quality Control Further information and registration at www.accuracy365.com

Accuracy 365 aims to monitor and compare, between laboratories, the results obtained in the Internal Quality Control processes of end-users of specific protein control materials manufactured by **TRIMERO Diagnostics**.

Main features:

- 100 % free and lifetime.
- Cloud application, Available 24 hours, 365 days per year.
- Intuitive interface and very simple to use.
- Calculation and representation of statistical data in real time.
- Anonymous participation, data are only accessible by its owner.

TRIMERO Diagnostics, SL
 c. València 558, 4t 2a - 08026 Barcelona (Spain)
 ☎ +34 93 244 86 79 - www.3diag.com

INSTRUCTIONS FOR USE

Reagents for professional use,
 for *In Vitro* use only in clinical laboratory (IVD)

3diag - sTfR - CONTROL

**Soluble Transferrin Receptor
 Controls (2 lev.)**

REF TD-42683

(Product included in **REF** TD-42680)

INTENDED USE

Internal controls, with assigned values, for the quantitative determination of Soluble Transferrin Receptor (sTfR), in human serum, by immunochemical methods.

PRINCIPLE OF THE METHOD

The specific antibodies (Ab) of the reagent, bound to polystyrene particles, when combined with the antigens (Ag) of the patient sample, form insoluble compounds causing a change in the absorbance and dispersion of the light, proportional to the antigen concentration, which can be quantified by turbidimetric (TIA) or nephelometric (NIA) method, by comparison with calibrators of known concentration.

CONTENTS - COMPOSITION - PREPARATION

- High Control: **CONTROL H sTfR**
REF TD-42683-H **CONT** 1 ml
- Low Control: **CONTROL L sTfR**
REF TD-42683-L **CONT** 1 ml

The controls are human serum solutions, delipidated, filtered by 0.2 µm.

As preservatives, the controls contain <0.1% (1 g/l) Sodium Azide.

The controls are ready for use and require no preparation.

Before each use it is convenient that the controls are homogenized, shaking them gently avoiding the formation of foam or bubbles

The values of the controls are lot dependent and are indicated in the table of values of their Instructions for Use.

WARNINGS - PRECAUTIONS

- Sodium Azide is toxic. Even if at the concentration present Sodium Azide is not harmful, take the necessary precautions to avoid accidental ingestion or contact with the eyes.
- Sodium Azide can react with lead or copper to give explosive compounds. For disposal it is recommended to rinse with plenty of running water to avoid accumulation in drains.
- Materials of human origin have been tested and found negative for the presence of HBsAg, HCV, and anti-HIV 1 and 2 antibodies.
- Since the absence of infectious agents can not be proven with absolute certainty, components containing materials of human or animal origin must be handled with caution, as potentially infectious, following the recommended safety standards for biological risk.

- Do not mix components belonging to different lot kits.
- Clinical diagnosis should not be based on the results of a single test, but should always integrate all relevant clinical and laboratory data.

STORAGE - SHELF LIFE

- Store refrigerated at +2...+8°C. Do not freeze, as the functionality of the controls may be altered.
- Properly stored and unopened, the controls are stable until the expiration date indicated on the label.
- Once opened, the shelf life of the controls is at least 4 weeks, provided that after each use they are stored immediately in the original containers, tightly capped and refrigerated at +2...+8°C. This information should be taken as a guideline given that, obviously, the shelf life depends on the particular environmental and use conditions, which may differ from those of the stability studies carried out.

MATERIALS NEEDED, NOT SUPPLIED

The controls are intended to be used in conjunction with the Reagents and Calibrators:

- 3diag - sTfR - 800** **REF** TD-42681
- 3diag - sTfR - TIA** **REF** TD-42691
- 3diag - sTfR - ADV** **REF** TD-42695
- 3diag - sTfR - ADV (5)** **REF** TD-42696
- 3diag - sTfR - CAL SET** **REF** TD-42682
- 3diag - sTfR - CAL** **REF** TD-42692

PROCEDURE

Follow the Instructions for Use of the analyzer used to program and calibrate an assay, with the general parameters recommended in the Instructions for Use of the Reagents.

In the analyzers of the *ADVIA® Chemistry Systems* series from *Siemens Healthcare* (*ADVIA®* and its associated brands are registered trademarks of *Siemens Healthcare*), to process the controls it is necessary to deactivate, for calibrators and controls, the analyzer's clot detection system. It may also be necessary to deactivate the clot detection system in other analyzers.

TRACEABILITY

Values are referred to the *Reference Reagent Recombinant Soluble Transferrin Receptor (rsTfR)* (NIBSC code: 07/202) of the WHO (*World Health Organization*), as monomers of free rsTfR⁽¹⁾.

(1) For more information visit the NIBSC website (*National Institute for Biological Standards and Control*): www.nibsc.org.

SYMBOLS

In addition to the harmonized symbols provided on the European Standard EN 980:2008, in the labels and instructions of use has been used the complementary symbology proposed⁽²⁾ by the *EDMA* (*European Diagnostic Manufacturers Association*), whose meaning is detailed below.

(2) EDMA Labelling Task Force: "*EDMA Symbols for IVD Reagents and Components - Revision, October 2009*".

CONTROL	Control
H	High
L	Low
sTfR	Soluble Transferrin Receptor
CONT	Contents

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
NOTE #1: Applications available for **Optilite®**, **Alinity c**, **Architect c**, **ADVIA®** series, **AU®** series, **cobas®** series, and other analyzers by request

IMMAGE® and **AU®** are registered trademarks of Beckman Coulter, Inc, Fullerton, CA.
BN™ and **ADVIA®** are registered trademarks of Siemens Healthcare Diagnostics Products GmbH, Marburg, Deutschland.
Optilite® is a registered trademark of The Binding Site Group Ltd, Birmingham, U.K.
Alinity and **Architect** related brand marks are registered trademarks of Abbott Laboratories, Abbott Park, Illinois, USA
cobas® and related brand marks are trademarks of Roche Diagnostics Ltd, Rottkreuz, Switzerland