
# C. difficile GDH FIA

STANDARD™ F C. difficile GDH FIA

SD BIOSENSOR

## KIT CONTENTS

## MATERIALS REQUIRED BUT NOT PROVIDED


- STANDARD F Analyzer

## SPECIMEN COLLECTION AND STORAGE

- Feces**
  - Fecal specimen should be collected in a clean and dry container and collected at any time of the day may be used.
  - Urine should be excluded from the fecal specimen.
  - Fresh fecal specimen should be used in this test. Do not use any transport media to store and transport the fecal specimen.
  - Fecal specimen may be stored in a refrigerator at 2-8°C / 35.8-46.4°F for up to 24 hours.
  - For prolonged storage, specimens may be frozen under -20°C / -4°F. It should be brought to room temperature prior to use.


## TEST PROCEDURE AND PREPARATION

- Preparation**
  - Allow kit components and collected specimen to room temperature at least 30 minutes before starting the test.
  - Carefully read instructions for the STANDARD F C. difficile GDH FIA.
  - Check the expiry date at the back of the foil pouch. Use another lot if expiry date has passed.
  - Open the foil pouch, and check the test device and desiccant inside the foil pouch.


## Collection of specimen

- Poke the 4 different areas of the fecal specimen using the sterile swab to collect a portion of feces (about 40-70mg). In case of liquid specimen, soak the sterile swab completely in the liquid fecal specimen.


The amount of fecal swab may affect the results. It is required to follow the swab amount of feces as shown in the picture above. Excessive fecal amount may induce a false positive result and slow migration.

- Take off the top of an extraction buffer tube and insert the swab into the tube. Swirl the swab at least ten times into the extraction buffer tube to dissolve the specimen.

- Remove the swab and dispose of the used swab in accordance with your biohazard waste disposal protocol. Tightly press down the filter cap onto the tube.

Do not use the extraction buffer tube without the filter cap. It affects the result value.

The composition of extraction buffer is same as STANDARD F C. difficile Toxin A/B FIAs.

## Analysis of specimen

### STANDARD TEST mode

- Take the test device out of foil pouch and place it on a flat and dry surface. Write patient information on the label of test device.

STANDARD F2400 analyzer "Workplace" -> "Run Test" -> Insert patient ID and / or operator ID on the analyzer  
STANDARD F200 analyzer "Standard Test" mode -> Insert patient ID and / or operator ID


## PERFORMANCE CHARACTERISTICS

### Analytical Performance


Source	Concentration
Bionote-U156	1.34 ng/ml

- Precision  
STANDARD F C. difficile GDH FIA totally meet the acceptance criteria of repeatability and reproducibility.


- Insert the test device to the test slot of the analyzer. When inserting the test device to the analyzer, the analyzer will read the barcode data, and check the test device is valid.

- Apply 3 drops of mixed specimen to the specimen well in the test device holding the prepared extraction buffer tube upside down.

- After applying the specimen, immediately press the "TEST START" button.


- The analyzer will automatically display the test result after 15 minutes.


### 'READ ONLY' mode

- Take the test device out of foil pouch and place it on a flat and dry surface. Write patient information on the label of test device.


- Prepare extracted specimens.


- Prepare test devices depending on the workload.

- Apply 3 drops of extracted specimen into test devices in sequence at about 30 seconds interval.


- Leave test devices for 15 minutes on a flat surface for incubation.

- Prepare a STANDARD F Analyzer and select the "Read Only" mode according to the analyzer's manual.

- Insert the test device to the analyzer.

- Do not put bubbles in the specimen well of the test device.
- Discard the used test kit in proper container, according to your healthcare professional's recommendation.

## INTERPRETATION OF TEST RESULT

Result	COI (Cutoff index) value	Interpretation
Positive	COI ≥ 1.0	Positive for <i>C. difficile</i> GDH
Negative	COI < 1.0	Negative for <i>C. difficile</i> GDH
Invalid	COI value is not displayed	Retest should be performed with a new test device and a new patient specimen

- Result must be interpreted together with STANDARD F *C. difficile* Toxin A/B FIA.
- False positive result can be caused by intake of antibiotics or other medicine.
- Retest should be performed with newly collected fecal specimen.

## EXPLANATION AND SUMMARY

### Introduction

*Clostridium difficile* is a gram-positive anaerobe. Infection with *C. difficile* causes severe diarrhea and can be fatal if not diagnosed and treated in a timely manner. *C. difficile* infection is induced in patients by long-term treatment with antibiotics and is commonly found in hospital environment. It is easily transmitted through contact with infected fecal matter. Accurate testing for toxins proved to be difficult due to assays having low sensitivity. Since all strains of *C. difficile* produce large amounts of glutamate dehydrogenase, testing for this antigen has proven to be a better screening tool due to its higher negative predictive value. *Clostridium difficile* Glutamate Dehydrogenase (GDH) is an enzyme produced in large quantities by all toxigenic and non-toxicogenic strains, making it an excellent marker for the organism. Diagnostic tests that detect GDH are highly sensitive. They have relatively low specificity for disease-causing *C. difficile* though, because GDH is also produced by non-disease-causing *C. difficile* strains and other *Clostridium* species.

### Intended use

STANDARD F C. difficile GDH FIA is a fluorescence immunoassay for the qualitative detection of *C. difficile* GDH from fecal specimens from persons. This is only for professional use. It provides only an initial screening test result. STANDARD F C. difficile GDH FIA should be used with the appropriate analyzer. STANDARD F Analyzers, manufactured by SD BIOSENSOR.

### Test Principle

STANDARD F C. difficile GDH FIA is based on the immunofluorescence technology with STANDARD F Analyzer manufactured by SD BIOSENSOR to measure the *C. difficile* GDH in human specimen. The specimen from human should be processed

for the preparation using the components of the STANDARD F C. difficile GDH FIA. After applying the specimen mixture to the test device, the complex will be formed on the membrane as the result of the antigen-antibody reaction. The intensity of the fluorescence light generated on the membrane is scanned by the STANDARD F Analyzer. STANDARD F Analyzers can analyze the *C. difficile* GDH of the clinical specimen based on a pre-programmed algorithms and display the test result on the screen.

## KIT STORAGE AND STABILITY

Store the kit at 2-30°C / 36-86°F, out of direct sunlight. Kit materials are stable until expiration date printed on the outer box. Do not freeze the kit.

## WARNINGS AND PRECAUTIONS

- Carefully follow the instructions and procedures described in this instructions before testing.
- STANDARD F C. difficile GDH FIA should be used with STANDARD F Analyzer.
- STANDARD F C. difficile GDH FIA should remain in its original sealed pouch until ready to use. Do not use if the pouch is damaged or the seal is broken.
- STANDARD F C. difficile GDH FIA is single use only. Do not re-use it.
- Do not use frozen specimen or artificial media.
- Place the analyzer on a flat surface when in use.
- Wash your hands in warm, soapy water. Rinse well and dry completely before testing.
- Discard the used test kit according to the proper method.
- Check the expiration date printed on the pouch or package.
- Use the STANDARD F C. difficile GDH FIA at 15-32°C / 59-60°F.
- All kit components must be at room temperature 30 minutes before running the assay.
- Do not write on the barcode or damage the barcode of the test device.

## QUALITY CONTROL

### Calibration

The calibration set test of STANDARD F Analyzers should be conducted according to the analyzer's manual.

### When to use calibration set

- Before using the analyzer for the first time
- When you drop the analyzer
- Whenever you do not agree with your result
- When you want to check the performance of an analyzer and test device

### How to use calibration set

- Select the "Calibration" menu.
- The specific calibration set is included with the analyzer.
- Insert the CAL-1 first, and then insert the CAL-2 for UV-LED testing and the CAL-3 for RGB-LED testing in order.

The STANDARD F analyzer automatically calibrate and identify the optical performance through measuring the membrane of the test device whenever the test is conducted in "Standard Test" mode. If "EEE" message displays on the screen, it means that the analyzer has a problem, so check with test device. Contact the SD BIOSENSOR local distributor if the "EEE" message still appears.

### External quality control

Quality control testing should be run to check the performance of STANDARD C. difficile GDH FIA and STANDARD F Analyzers. STANDARD C. difficile GDH FIA manufactured by SD BIOSENSOR should be used for quality control testing. Control test should be conducted in accordance with the instructions of STANDARD C. difficile GDH Control.

Control test should be run:

- once for each new lot.
- once for each untrained operator.
- as required by test procedures in instructions for use of STANDARD C. difficile GDH Control and in accordance with local, state and federal regulations or accreditation requirements.

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## BIBLIOGRAFIA

1. Wren, M.W.D, et al. "Laboratory diagnosis of infection. An evaluation of tests for faecal toxin, glutamate dehydrogenase, lactoferrin and toxigenic culture in the diagnostic laboratory". British Journal of Biomedical Science, 66 (1), 2009.
2. Vaishnavi, Ch., "Clinical spectrum & pathogenesis of difficile associated diseases". Indian J. Med. Res. 131, April 2010, pp 487-499.
3. Sharp SE, Ruden LO, Pohj J, et al. Evaluation of the Cd Diff Quik Check Complete Assay, a new glutamate dehydrogenase and A/B toxin combination lateral flow assay for use in rapid, simple diagnosis of Clostridium difficile disease. J Clin Microbiol. 2010;48:2082-2086.
4. Surawicz CM, Brandt LJ, Binion DG, et al. Guidelines for diagnosis, treatment, and prevention of Clostridium difficile infections. Am J Gastroenterol. 2013;108:478-498.
5. Reller ME, Alcabasa RC, Lema CA, et al. Comparison of two rapid assays for Clostridium difficile common antigen and a C difficile toxin A/B assay with the cell culture neutralization assay. Am J Clin Pathol. 2010;133:107-109.

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- STANDARD F-Analysegerät

## SYMBOL

REF	Reference number	Caution	Use by	LOT	Batch code	Consult Instructions for Use	Do not re-use
IVD	In vitro Diagnostics	Note	Manufacturer		Date of manufacture	Contains Sufficient for <math>n</math>-Tests	Keep away from sunlight
	Indicate that you should keep the product dry	To indicate the temperature limitations in which the transport package has to be kept and handled.			This product fulfills the requirements of the European Directive 98/79/EC.		Do not use if packaging is damaged

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