


# EN

## STANDARD F COVID-19 Ag FIA

STANDARD™ F COVID-19 Ag FIA

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

### KIT CONTENTS




### MATERIALS REQUIRED BUT NOT PROVIDED

• STANDARD F Analyzer

• Timer

### SPECIMEN COLLECTION AND STORAGE

#### Nasopharyngeal swab



- To collect a nasopharyngeal swab specimen, insert a sterile swab into the nostril that presents the most secretions.
- While rotating the swab, insert swab approximately one inch (about 2 cm) into nostril until resistance is met at turbinates.
- Gently rotate the swab 3 to 4 times against nasal/pharyngeal wall, and leave swab in place for several seconds to absorb secretions.
- Remove the sterile swab from the nostril carefully.
- Specimens should be as soon as possible after collection.
- Specimen in the extraction buffer tube may be stored at room temperature for up to 1 hours or at -20°C/-36°F for up to 40 days.

#### Nasal swab



- Tilt patient's head back slightly.
- While rotating the swab, insert swab approximately one inch (about 2 cm) into nostril until resistance is met at turbinates.
- Slowly rotate the swab 3 to 4 times against the nasal wall 5 times for a minimum of 15 seconds.
- Remove the swab from the nostril carefully.
- Specimens should be as soon as possible after collection.
- Specimen in the extraction buffer tube may be stored at room temperature for up to 1 hours or at -20°C/-36°F for up to 40 days.

#### Viral transport medium

- The following viral transport media listed in Table 1 were validated by SD BIOSENSOR's R&D Department and determined to be compatible with the STANDARD F COVID-19 Ag FIA. The viral transport media are divided into two different groups by manufacturing brand. Find a brand of viral transport medium you use, and select a correct group button while sample selecting phase.

Table 1. Recommended Viral Transport Media	
<b>Group</b>	<b>Viral Transport Medium</b>
VTM - Group 1	COPAN UTM™ STANDARD™ Transport Medium BD Universal Viral Transport
VTM - Group 2	


- Specimens in transport media should be transported directly to the laboratory, and be preferably processed immediately. If immediate delivery or processing is delayed, the specimen should be stored at -20°C and processed within 8 hours. If delivery or processing exceed temperature conditions above, specimens should be transported in dry ice and once in laboratory frozen at -20°C or colder.
- Only viral transport medium verified by SD BIOSENSOR can be used as specimen for STANDARD F COVID-19 Ag FIA. Please check the brand of viral transport medium and select a correct group button while sample selecting phase.

### PREPARATION AND TEST PROCEDURE

#### ■ Preparation

- Allow test device and selected specimen to room temperature (15-30/59-86°F) prior to testing.
  - Check the expiry date at the back of the foil pouch of the STANDARD F COVID-19 Ag FIA.
  - Check the condition of the test device and desiccant before use.
- The violet-line on the membrane of unused test device will disappear after use.  
Do not write on the barcode or damage the barcode of the test device.

- Extraction of specimen [Nasopharyngeal and Nasal swab]



- Place the swab specimen into the extraction buffer tube. While squeezing the buffer tube, stir the swab more than 5 times.
- Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.
- Press the nozzle cap tightly onto the tube.

#### [Specimen in viral transport medium]



- Collect 350 µl of viral transport medium specimen using a micro pipette.
- Add 350 µl of viral transport medium specimen into the extraction buffer tube and press the nozzle cap tightly onto the tube.

#### PERFORMANCE CHARACTERISTICS

- Clinical evaluation  
To evaluate the clinical sensitivity and specificity of STANDARD F COVID-19 Ag FIA, specimens were collected from those who had positive and negative results of COVID-19 infection through RT-PCR method. 160 specimens including specimens with RT-PCR Ct value >30 and 100 negative specimens were tested at the hospital laboratory.

#### 1) Clinical Sensitivity (C>30) and specificity : Nasal swab

	RT-PCR	Positive	Negative	Total
STANDARD F COVID-19 Ag FIA	Positive	49	0	49
	Negative	3	100	103
Total	52	100	152	

#### Clinical sensitivity

94.23% (95% CI: 84.05% - 98.79%)

The positive predictive value(PPV) was 100.00%

The negative predictive value(NPV) was 97.09% (95% CI: 91.74% - 99.01%)

The sensitivity including 8 patients with RT-PCR Ct value >30 was 66.67%

#### 2) Clinical Sensitivity (C>30) and specificity : Nasopharyngeal swab

	RT-PCR	Positive	Negative	Total
STANDARD F COVID-19 Ag FIA	Positive	49	0	49
	Negative	3	100	103
Total	52	100	152	

#### Clinical sensitivity

94.23% (95% CI: 84.05% - 98.79%)

The positive predictive value(PPV) was 100.00%

The negative predictive value(NPV) was 97.09% (95% CI: 91.74% - 99.01%)

The sensitivity including 8 patients with RT-PCR Ct value >30 was 66.67%

#### RT-PCR

	Positive	Negative	Total	
STANDARD F COVID-19 Ag FIA	Positive	49	0	49
	Negative	3	100	103
Total	52	100	152	

#### Clinical specificity

100.00% (95% CI: 96.38% - 100.00%)

#### ■ Clinical performance

The positive predictive value(PPV) was 100.00%

The negative predictive value(NPV) was 97.09% (95% CI: 91.74% - 99.01%)

The sensitivity including 8 patients with RT-PCR Ct value >30 was 66.67%

#### ■ Cross-reactivity

No cross-reactivity was observed for the following microorganisms at the indicated concentrations. All microorganisms were spiked into negative control matrix for testing.

#### ■ Intended use

STANDARD F COVID-19 Ag FIA is a fluorescent immunoassay for the qualitative detection of the specific nucleic acid sequence of SARS-CoV-2 in nasopharyngeal and nasal swab specimens.

STANDARD F COVID-19 Ag FIA should be used with the STANDARD F Analyzers manufactured by SD BIOSENSOR. This test is for *in vitro* professional diagnostic use and is intended as a tool for early diagnosis of SARS-CoV-2 infection in patient with clinical symptoms with SARS-CoV-2 infection. It provides only an initial screening result. For specific alternative diagnostic methods should be performed in order to obtain the confirmation of SARS-CoV-2 infection.

#### ■ External quality control

1) STANDARD COVID-19 Ag Control may be used to validate that reagents and assay procedure perform properly.

2) STANDARD COVID-19 Ag Control is available to be purchased separately.

3) It is recommended that positive and negative controls to run for each untrained operator.

- once for each new lot.

- as required by test procedures in this instructions and in accordance with local.

#### INTERPRETATION OF TEST RESULT

Result	COI (Cut-off index) value	SARS-CoV-2 Ag
Positive	COI > 1.0	Positive for SARS-CoV-2 Ag
Negative	COI < 1.0	Negative for SARS-CoV-2 Ag
Invalid	COI value is not displayed	Retest should be performed with a new test device and a new patient's specimen.

The test result of a specimen is given either as Positive (+) or Negative (-). The COI is a numerical representation of the measured fluorescence signal.

Positive results should be considered in conjunction with the clinical history and other data available to the physician.

#### INTERPRETATION OF THE TEST

Result	COI (Cut-off index) value	SARS-CoV-2 Ag
Positive	COI > 1.0	Positive for SARS-CoV-2 Ag
Negative	COI < 1.0	Negative for SARS-CoV-2 Ag
Invalid	COI value is not displayed	Retest should be performed with a new test device and a new patient's specimen.

The test result of a specimen is given either as Positive (+) or Negative (-). The COI is a numerical representation of the measured fluorescence signal.

Positive results should be considered in conjunction with the clinical history and other data available to the physician.

#### INTERPRETATION DE LA PRÉPARATION Y PRUEBA

##### ■ Preparación

1. Take the test device out of the foil pouch and place it on a flat and dry surface. Write patient information on the label of the test device.

2. Prepare extracted specimen.

3. Prepare test devices depending on the workload.

4. Apply 4 drops of extracted specimen to the specimen well of the test device.

5. After applying the specimen, immediately press TEST START button.

6. The analyzer will automatically display the test after 15 minutes.

#### INTERPRETATION DE LOS RESULTADOS

Resultado	Valor COI (índice de señal)	SARS-CoV-2 Ag
Positivo	COI > 1.0	Positivo para SARS-CoV-2 Ag
Negativo	COI < 1.0	Negativo para SARS-CoV-2 Ag
Inválido	No se muestra valor COI	Se debe realizar una prueba con un nuevo dispositivo de prueba y una nueva muestra del paciente.

El resultado de la prueba de una muestra se da como Positivo (+) o Negativo (-).

(- Neg (-) con un valor COI (del inglés Cut off Index). El COI es una representación numérica de señal de fluorescencia medida.

Los resultados positivos deben ser considerados en conjunto con el historial clínico y demás información de apoyo al diagnóstico de la prueba.

#### INTERPRETACIÓN DE LOS RESULTADOS DE TEST

Résultat	Valeur COI (Indice de seuil)	SARS-CoV-2 Ag
Positif	COI > 1.0	Positif sur SARS-CoV-2 Ag
Négatif	COI < 1.0	Négatif pour SARS-CoV-2 Ag
Inválido	La valeur COI n'est pas affichée	Il faut réaliser le test avec un nouveau dispositif de test et un nouvel échantillon prélevé sur un patient.

Le résultat d'un échantillon est soit positif (+) soit négatif (-) avec une valeur COI (ou indice de seuil).

