**EN** FEP F-CRP

**STANDARD F CRP**  
STANDARD™ F CRP  
PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

**SD BIOSENSOR**

**KIT CONTENTS**

Test device (individually in a foil pouch with desiccant and Spoit™)  
Extraction buffer  
Instructions for use

**MATERIALS REQUIRED BUT NOT PROVIDED**

- STANDARD F Analyzer

**SPECIMEN COLLECTION AND STORAGE**

- Whole blood**  
(Capillary whole blood)
- Capillary whole blood should be collected aseptically by fingertip.
  - Clean the area to be lanced with an alcohol swab.
  - Squeeze the end of the fingertip and pierce with a sterile lancet.
  - Collect the accurate volume of capillary whole blood using the Spoit (Red) in kit for the testing.
  - The capillary whole blood must be tested immediately after collection.

- Venous whole blood**
- Collect the venous whole blood into the commercially available anticoagulant tube such as heparin or EDTA by venipuncture.
  - It is recommended that collected venous whole blood specimens are used immediately. If venous whole blood in an anticoagulant tube is stored in a refrigerator at 2-8°C/36-46°F, the specimen can be used for testing within 8 hours after collection.
  - Do not use hemolyzed blood specimens.


- Serum**
- Collect the whole blood into the commercially available plain tube, NOT containing anticoagulants such as heparin or EDTA by venipuncture and leave to settle for 30 minutes for blood coagulation and then centrifuge blood to get serum specimen of supernatant.
  - Serum in the plain tube may be stored at room temperature (15-30°C/59-86°F) for up to 8 hours and at 2-8°C/36-46°F for up to 3 days prior to testing.
  - For over 3 days storage, specimens may be frozen under -40°C/-40°F for up to 3 months.
  - It should be brought to room temperature prior to use.

- Plasma**
- Collect the venous whole blood into the commercially available anticoagulant tube such as heparin or EDTA by venipuncture and centrifuge blood to get plasma specimen.
  - Plasma in an anticoagulant tube may be stored at room temperature (15-30°C/59-86°F) for up to 8 hours and at 2-8°C/36-46°F for up to 3 days prior to testing.
  - For over 3 days storage, specimens may be frozen under -40°C/-40°F for up to 3 months.
  - It should be brought to room temperature prior to use.

- Anticoagulants such as heparin or EDTA do not affect the test result.
- As known relevant interference, hemolytic specimens, rheumatoid factors-contained specimens and lipemic, icteric specimens can lead to impaired the test results.

**PREPARATION AND TEST PROCEDURE**

- Preparation**
- Allow kit components and collected specimen to room temperature (15-30°C / 59-86°F) at least 30 minutes before starting the test.
  - Carefully read instructions for the STANDARD F CRP.
  - Check the expiry date at the back of the foil pouch. Use another lot if expiry date has passed.
  - Open the foil pouch, and check a test device, desiccant and a Spoit containing a CRP latex tablet.

Do not write on the barcode or damage the barcode of the test device.

**Analysis of specimen**

STANDARD TEST™ mode	STANDARD F200 and F2400 analyzer
---------------------	----------------------------------

- Take the test device out of the foil pouch and place it on a flat and dry surface. Write patient information on the label of test device.

STANDARD F2400 analyzer	Workplace → Run Test → Scan or type patient ID and/or operator ID
STANDARD F200 analyzer	Standard Test™ mode → Insert patient ID and / or operator ID on the analyzer

- Insert the test device to the test slot of the analyzer. When inserting the test device to the analyzer, the analyzer read the barcode data, and check the test device is valid.

- Collect S-J of blood specimen with the Spoit. The specimen will be collected automatically by the capillary action when the tip of the Spoit touches the droplet of the blood.

**QUALITY CONTROL**

- Calibration**
- Before using the analyzer for the first time.
  - When you drop the analyzer.

**PERFORMANCE CHARACTERISTICS**


**ANALYTICAL PERFORMANCE**

**1. Analytical Sensitivity – LoB, LoD and LoQ**

The analytical sensitivity (Limit of Blank (LoB), Limit of Detection (LoD) and Limit of Quantitation (LoQ)) of STANDARD™ F CRP is shown below;

STANDARD™ F CRP	C-reactive protein (mg/L)		
	LoB	LoD	LoQ
Serum	0.42	0.93	1.0
Whole Blood	0.42	0.97	1.0

- Put the tip of the Spoit into an extraction buffer. Mix the collected specimen, CRP latex and extraction buffer by carefully pressing and releasing the rubber at the top of the Spoit for 6-8 times. Mix carefully to avoid bubble forming.
- Collect all the specimen mixture with the Spoit from the tube and apply to the specimen well of the test device.
- After applying the specimen, immediately press the TEST START button.
- The analyzer will automatically display the test result after 3 minutes.

- Do not put bubbles in the specimen well of the test device.
- Discard the used test kit in proper container according to your healthcare professional's recommendation.

**INTERPRETATION OF TEST RESULTS**

**Introduction**

STANDARD F Analyzers used with STANDARD F CRP reads CRP concentration between 1.0-150mg/L for capillary or venous whole blood specimen and 1.0-130mg/L for plasma or serum specimen. If the result is below 1.0mg/L for whole blood, serum or plasma it will be reported as < 1.0mg/L. If the result is above 150mg/L for whole blood and 130mg/L for serum of plasma, it will be reported as > 150mg/L and > 130mg/L each.

**\*Reference value: ≥ 10mg/L**  
10mg/L is generally considered as a clinical cut-off value for inflammatory disease.\*

**EXPLANATION AND SUMMARY**

**Intended use**

STANDARD F CRP is an *in vitro* diagnostic use to measure the CRP in the human serum, plasma and whole blood specimen. The measurement of CRP provides information for the detection and evaluation of infection, tissue injury, inflammatory disorders and associated diseases. STANDARD F CRP should be used with the appropriate analyzer, STANDARD F Analyzers, manufactured by SD BIOSENSOR.

**Test Principle**

STANDARD F CRP is based on the reflectometry and immunosay technology with STANDARD F Analyzer manufactured by SD BIOSENSOR. The specimen from human should be processed for the preparation using the components of the STANDARD F CRP. CRP in human specimen is react with monoclonal anti-CRP, latex in the tablet of the Spoit (Red) and makes the complexes. When applying the specimen mixture to the test device, the complexes move along the membrane and are captured by the monoclonal anti-CRP in the test line on the membrane as the result of the antigen-antibody reaction. The intensity of the reflected light is scanned and converted into an electric signal which is proportional to the intensity of the light produced on the membrane. STANDARD F Analyzers can analyze the CRP concentration of the clinical specimen based on preprogrammed algorithms and display the test result on the screen.

**KIT STORAGE AND STABILITY**

Store the kit at 2-30°C / 36-86°F out of direct sunlight. Kit materials are stable until expiration date printed on the outer box. Do not freeze the kit.

**WARNINGS AND PRECAUTIONS**

- STANDARD F CRP is for *in vitro* diagnostics use only.
- Carefully follow the instructions and procedures described in this instructions before testing.
- Do not use extraction buffer of another lot.
- Use the STANDARD F CRP at 15-32°C / 59-90°F and 10-90%RH.
- All kit components and collected specimens are must be at room temperature at least 30 minutes before printing the assay.
- Check the expiration date printed at the pouch or package. Do not use after expiration date.
- Do not use the test kit if the pouch is damaged or the seal is broken.
- Only single use. Do not re-use it.
- Observe established precautions against microbiological hazards throughout testing procedures.
- Dispose of all specimens and materials used to perform the test as bio-hazard waste. Laboratory chemical and biohazard wastes must be handled and discarded F Analyzers can analyze the CRP concentration of the clinical specimens based on preprogrammed algorithms and display the test result on the screen.
- Descartar in foil pouch is to absorb moisture and keep humidity from affecting products. If the moisture indicating desiccant beads change from yellow to green, the test device in the pouch should be discarded.
- Wash your hands in warm and soapy water. Rinse well and dry completely before testing.
- Do not use hemolyzed specimens.
- Do not use frozen whole blood specimen or any artificial materials.
- Once a test device has been successfully scanned by analyzer, do not attempt to scan the test device again in the same analyzer.
- Improper specimen collection, handling or transportation may yield inaccurate results.
- Do not write on the barcode or damage the barcode of the test device.

**LIMITATION OF TEST**

- The test should be used for the detection of CRP levels in human whole blood, serum and plasma specimens.
- Failure to follow the test procedure and interpretation of test results may adversely affect test performance and/or produce invalid results.
- Invalid results may occur if a poor-quality specimen is obtained.
- The test result must always be evaluated with other clinical data available to the physician.

**QUALITY CONTROL**

- Calibration**
- Before using the analyzer for the first time.
  - When you drop the analyzer.

**PERFORMANCE CHARACTERISTICS**

**2. Precision**

The following precision results of STANDARD™ F CRP meet the acceptance criteria of both repeatability and reproducibility (CV ≤15%).

CRP (mg/mL)	Within Run									
	Level 1			Level 2			Level 3			
	Site 1	Site 2	Site 3	Site 1	Site 2	Site 3	Site 1	Site 2	Site 3	
Serum	AVG.	11.5	11.6	11.7	48.4	48.4	48.9	102	104	100
	CV (%)	11.8	12.1	11.9	49.5	50.3	51.3	105	108	103
	DIF (%)	2.3%	4.5%	1.8%	2.4%	3.9%	4.9%	2.6%	4.0%	3.2%
Plasma	AVG.	10.5	10.4	10.3	49.8	50.4	49.5	105	107	103
	CV (%)	10.9	10.8	10.6	50.8	52.0	51.6	109.5	112.0	106.2
	DIF (%)	8.4%	7.0%	7.0%	9.0%	7.6%	5.4%	8.4%	6.9%	6.5%

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