**STANDARD F
iFOB FIA**

STANDARD™ F iFOB FIA

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

SD BIOSENSOR

KIT CONTENTS**MATERIALS REQUIRED BUT NOT PROVIDED**


• STANDARD F Analyzer

SPECIMEN COLLECTION AND PREPARATION**Feces**

1. Fecal specimen should be collected in a clean and dry container and collected at room temperature.
2. Urine should be excluded from the fecal specimen.
3. Fresh fecal specimen should be used in this test. Do not use any transport media to store and transport the fecal specimen.
4. Fecal specimen may be stored at room temperature for up to 24 hours or at 2-8°C/36-46°F for up to 15 days prior to testing.

TEST PROCEDURE**Preparation**

1. Allow kit components and collected specimen to room temperature at least 30 minutes before starting the test.
2. Carefully read instructions for the STANDARD F iFOB FIA.
3. Check the expiry date at the back of the foil pouch. Use another lot if expiry date has passed.
4. Open the foil pouch and check the test device in the foil pouch.

Before Use After Use
Do not write on the barcode or damage the barcode of the test device.

Specimen collection

1. Open the green cap of the Sample collection tube.
2. Stab 5 to 6 sites of the fecal specimen using the stick (Do not scoop feces).

3. Close the green cap tightly, and then shake the tube several times vigorously.

- Specimen mixture (assay diluent + feces) is able to leak if the white cap is opened. Please do not open the white cap while collecting the specimen.
- Collected sample tube can be stored in room temperature for 24 hours. If you need to store it more than 24 hours, please store it in the refrigerator (2-8°C/36-46°F) for up to 15 days.

Analysis of specimen**'STANDARD TEST' mode****STANDARD F200 and F2400 analyzer**

1. Take the test device out of the foil pouch and place it on a flat and dry surface. Write patient information on the label of the test device.

STANDARD F2400 analyzer	'Workplace' → 'Run Test' → Scan or type patient ID and/or operator ID
STANDARD F200 analyzer	'Standard Test' mode → Insert patient ID and/or operator ID on the analyzer

2. Insert the test device to the test slot on the analyzer. After inserting the test device to the analyzer, the analyzer read the barcode data, and check the test device is valid.

3. Shake the Sample collection tube several times, and open the white cap.

4. Dispense 3 drops of specimen mixture into the specimen well of the test device.

PERFORMANCE CHARACTERISTICS**Analytical performance****1. Accuracy (Method Comparison)**

The test results of the fecal samples comparing the STANDARD™ F iFOB FIA with an automated immunochemical method are shown below;

Regression Analysis

Slope: 1.005

Y-Intercept: 3.582

R: 0.976

R²: 0.952**Acceptable System Accuracy**

Below -35%: 2 / 210 (1%)

Within ±35%: 203 / 210 (96.6%)

Over +35%: 5 / 210 (2.4%)

5. After applying the specimen, immediately press the 'TEST START' button.

6. The analyzer will automatically display the test result after 5 minutes.

**INTERPRETATION OF TEST RESULTS****Measuring Range**

25 - 1,000ng/mL (5 - 200 µg Hb/g feces)

- If below 25 ng/mL, '25-' message will be displayed.

- If above 1,000 ng/mL, '1,000+' message will be displayed.

Out of measuring range

< 100 ng/mL (20µg Hb/g feces)

*It is possible to adjust cut-off concentration according to study settings, national guidelines, available local resources, as well as age and gender of subjects. Most studies recommend a cut-off concentration of 75-1000ng/mL (15-200 µg Hb/g feces).

Clinicians should use the test results in conjunction with the patient's other diagnostic findings and clinical signs, and interpret the concrete values in the context of the patient's clinical situation.

EXPLANATION AND SUMMARY**Introduction**

Fecal occult blood (FOB) test is conducted to detect the occult blood in feces namely a quantity of hemoglobin unidentified with the naked eye. Examination of FOB is important in the case of diagnosis and treatment of many diseases including gastric ulceration, digestive tract cancer, tuberculosis, dysentery, and epidemic typhus. Especially, FOB test is a strong useful guide for medical treatment since a positive FOB result appears in about 50-70% of gastric ulcer patients and more than 80% of colorectal cancer patients. However, even though the result of FOB test is negative, the probability of relevant diseases cannot exclude so that it is recommended to conduct the repeat exams through the time and the diagnosis of diseases should be concluded with other laboratory test and the opinion of medical team. STANDARD F iFOB FIA is a simple and inexpensive non-invasive test gives an advantage to be able to screen the relative diseases.

Intended use

STANDARD F iFOB FIA is a fecal immunochemical test (FIT) for quantitative detection of hemoglobin (Hb) in fecal specimen. Quantitative FOB test can be used for detecting excess bleeding from the lower gastro-intestinal channel. (Inflammatory bowel disease, colorectal cancer screening programs). STANDARD F iFOB FIA provides a reliable and convenient tool for fecal occult blood quantification in both laboratory and point-of-care.

Test Principle

STANDARD F iFOB FIA is based on the immunofluorescence technology with STANDARD F Analyzer manufactured by SD BIOSENSOR to measure the hemoglobin concentration in fecal specimen. The test device consists of the test tube and the control tube. After the preparation using the components of the STANDARD F iFOB FIA, After applying the specimen mixture to the test device, the complex will be formed on the membrane as the result of the antigen-antibody reaction. The intensity of fluorescence light is measured on the membrane. STANDARD F Analyzers can analyze the hemoglobin concentration of the clinical specimen based on a pre-programmed algorithms and display the test result on the screen.

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Kit Storage and Stability

For 18 months when stored at 2-30°C/36-86°F and 10-90% RH.

WARNINGS AND PRECAUTIONS**1. Standard F iFOB FIA is for in vitro diagnostics use only.**

2. Carefully follow the instructions and procedures described in this Instructions for use before testing.

3. STANDARD F iFOB FIA should be used with STANDARD F Analyzer.

4. STANDARD F iFOB FIA should remain in its original sealed pouch until ready to use.

5. STANDARD F iFOB FIA is single use only. Do not re-use it.

6. Do not use any artificial materials as specimen.

7. Place the analyzer on a flat surface when in use.

8. Wash your hands in warm, soapy water. Rinse well and dry completely before testing.

9. Discard the used test kit according to the proper method.

10. Desiccant in foil pack is to absorb moisture and keep humidity from affecting products. If the moisture indicating desiccant beads change from yellow to green, the test device in the pouch should be discarded.

11. Check the expiration date printed at the pouch or package.

12. Use a volume (2mL) of sample collection tube with assay diluent.

13. Use the STANDARD F iFOB FIA at room temperature.

14. All kit components must be at room temperature 30 minutes before running the assay.

15. Do not write on the barcode or damage the barcode of the test device.

LIMITATION OF TEST

1. The test should be used for the detection of fecal occult blood levels in human fecal specimen.

2. Failure to follow the test procedure and interpretation of test results may adversely affect test performance and/or produce invalid results.

3. Invalid results may occur if a poor-quality specimen is obtained.

4. The test result must always be evaluated with other clinical data available to the physician.

Quality Control**Calibration**

The calibration set of STANDARD F Analyzer should be conducted according to the analyzer's manual.

[When to use calibration set]

1. Before using the analyzer for the first time.
2. When you drop the analyzer.
3. Whenever you do not agree with the final result.
4. When you want to check the performance of the analyzer and test device.

[How to use calibration set]

Calibration set is a required function that ensures optimal performance by checking the internal analyzer optics and functions.

1. Select the 'Calibration' in main menu.

2. The specific calibration set is included with the analyzer.

3. Insert the CAL-1 for white calibration, CAL-2 for UV LED calibration, and CAL-3 for RGB LED calibration in sequence.

STANDARD F Analyzer automatically calibrate and the optical performance of the test device whenever the test is conducted in 'Standard Test' mode. If the 'EEE' message displays on the screen, it means that the analyzer has a problem, so check with test device. Contact the SD BIOSENSOR local distributor if the 'EEE' message still appears.**Internal quality control**

1. The internal procedural control zone is on the membrane of the test device. STANDARD F Analyzers can automatically identify the internal procedural control zone and decide the result is valid or invalid.

2. The invalid result denotes that the signal is not within the pre-set range. If the screen of STANDARD F analyzer shows 'Invalid', turn off and turn on the analyzer again and re-test with a new test device.

External quality control
External quality control should be run to check the performance of STANDARD F iFOB FIA and STANDARD F Analyzer. STANDARD F iFOB FIA can be used for the external quality control. Control test should be conducted in accordance with the Instructions for use of STANDARD F iFOB Control. External quality control test should be run:

- once a month
- once for each untrained operator.
- as required by test procedures in the Instructions for use of STANDARD F iFOB Control and in accordance with local, state and federal regulations or accreditation requirements.

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