PIPETMAN® F

ENGLISH

GILSON®
SOLUTIONS AT WORK FOR YOU

PIPETMAN® F

CONTENTS	page
1 - INTRODUCTION	2
2 - SPECIFICATIONS	3
3 - OPERATION	4
4 - USER ADJUSTMENT	5
5 - PRE-RINSING	6
6 - GILSON DIAMOND® TIPS	7
7 - RECOMMENDATIONS	9
8 - TROUBLESHOOTING	10
9 - REPLACEMENT PARTS	11

1 - INTRODUCTION

Pipetman® F fixed volume pipettes are available in 13 different models. Together, they cover a range of volumes from 2 to 1000 μ L.

The pipetting mechanism incorporates a highly polished stainless steel piston, a seal, and an O-ring. The mechanism requires no routine lubrication. The high impact fluorocarbon moulded handgrip (Fig. 1D) has a low thermal coefficient.

The Pipetman F uses disposable polypropylene tips (Fig. 1B) to avoid cross contamination between samples. To protect the user from contamination by the tips, the Pipetman F is equipped with a built-in tip-ejector (Fig. 1A).

The tip-holder (Fig. 1C), the connecting nut and the tip-ejector are autoclavable for 20 minutes


Figure 1

at 121°C and 0.1 MPa pressure.

The Pipetman F can be easily re-adjusted in the laboratory, with compensation for variations in fluid viscosity and density. A special calibration key is provided.

The nominal volume of each model is shown on the small label on the push-button (Fig. 1E).

2 - SPECIFICATIONS

Pipetman F is a high quality pipette that offers excellent accuracy and precision. The figures given in the 'Gilson Maximum Permissible Errors' table were obtained using 'Gilson Diamond Tips'. These figures are only guaranteed using Genuine Gilson Diamond Tips. Pipetman F2 to F200 are tested using D200 Diamond Tips and F200 to F1000 are tested using D1000.

Gilson Maximum Permissible Errors

Reference	Model	Accuracy (systematic error)		Precision (random error)	
		Absolute μL	Relative %	Absolute S. D. μL	Relative %
F123770	F 2	± 0.10	± 5.00	≤ 0.03	≤ 1.50
F123771	F 5	± 0.10	± 2.00	≤ 0.04	≤ 0.80
F123772	F 10	± 0.10	± 1.00	≤ 0.05	≤ 0.50
F123604	F 20	± 0.20	± 1.00	≤ 0.06	≤ 0.30
F123775	F 25	± 0.25	± 1.00	≤ 0.07	≤ 0.30
F123778	F 50	± 0.40	± 0.80	≤ 0.15	≤ 0.30
F123784	F 100	± 0.80	± 0.80	≤ 0.25	≤ 0.25
F123605	F 200	± 1.60	± 0.80	≤ 0.30	≤ 0.15
F123787	F 250	± 3.00	± 1.20	≤ 0.75	≤ 0.30
F123788	F 300	± 3.50	± 1.17	≤ 0.75	≤ 0.25
F123789	F 400	± 3.60	± 0.90	≤ 0.80	≤ 0.20
F123790	F 500	± 4.00	± 0.80	≤ 1.00	≤ 0.20
F123606	F 1000	± 8.00	± 0.80	≤ 1.30	≤ 0.13

 *The data given in the table conform to the ISO 8655-2 Standard.*

Each pipette is inspected and validated by qualified technicians according to the Gilson Quality System.

PIPETMAN® F

Gilson declares that its manufactured pipettes comply with the requirements of the ISO 8655 Standard, by type testing. The adjustment is carried out under strictly defined and monitored conditions (ISO 8655-6):

- Basis of adjustment, Ex.
- Reference temperature, 20 °C
- Relative humidity, 50 %
- Barometric pressure, 101 kPa
- Use of distilled water grade 3 (ISO 3696)
- Ten measurements for each test volume, which are Nominal Volume, 50 % of Nominal Volume, and the minimum or 10 % of Nominal Volume.

3 - OPERATION

Place a tip on the tip-holder. Press the tip on firmly using a slight twisting motion to ensure a positive, airtight seal. (See Chapter 6 for the appropriate tip.)

 *Always fit a tip to your Pipetman before aspirating any liquid.*

Aspirating

- Press the push-button to the first positive stop (Fig. 2A).
- Holding the pipette vertically, immerse the tip into the liquid (depth immersion 2-4 mm).
- Release the push-button **slowly** and **smoothly** to aspirate the liquid (Fig. 2B) ;
- Wait one second and then withdraw the tip from the liquid. Wipe any droplets away from the outside of the tip using a medical wipe.

 *Avoid touching the orifice of the tip.*

Dispensing

- Place the end of the tip against the inside


Figure 2

wall of the vessel at an angle of 10 to 40 degrees. Press the push-button **smoothly** to the first stop (Fig. 2C). Wait one second. Press the push-button to the second stop to expel any remaining liquid (Fig. 2D).

- Keeping the push-button pressed to the end, remove the pipette by drawing the tip along the inside surface of the vessel.
- Release the push-button.
- Eject the tip by pressing the tip-ejector button. It is only necessary to change the tip if a different liquid is being sampled.

4 - USER ADJUSTMENTS

Pipetman F is factory calibrated using distilled water and very high precision balances. The nominal value of the Pipetman F may be slightly adjusted to compensate for liquids of different density or viscosity.

Adjustment of the factory calibration must only be performed using a special calibration key (Figures 3 and 5G).

This adjustment to accommodate for density or viscosity is limited to one full turn of the calibration key in either direction, which equals:

- ± 1.0 μL (models F2 to F20)
- ± 4.0 μL (models F25 to F100)
- ± 10.0 μL (models F120 to F200)
- ± 40.0 μL (models F250 to F1000).

 *The display is graduated in tenths of the key revolution (letters A, B,... J) and is read from left to right. Each small graduation mark represents 0.2 of the distance between each letter.*

To adjust the volume setting to compensate for a specific density or viscosity, slide the key over the push-button until notches fit into the grooves in the calibration sleeve (Figure 3), then turn the key:

PIPETMAN® F

- clockwise to decrease the volume;
- counterclockwise to increase the volume; go above the desired value, then decrease the volume to reach the required setting.

Example: When using a Pipetman F100 to aspirate a particularly viscous solution, you may determine gravimetrically that the volume delivered is 98.8 μL , and the display reads 'E 0.2'. Therefore, you wish to increase the volume dispensed by 1.2 μL . As the interval between each letter on the calibration display represents 0.4 μL for the F100 (1/10 of 4 μL), the calibration tool must be turned 3 units until the display reads 'B 0.2'.


Figure 3

It is possible to alter the volume by more than one full turn in the counterclockwise direction. In this case, remember to readjust the same number of turns when it is necessary to return to the original setting.

Check the new volume gravimetrically. If the volumes delivered are still not sufficiently close to the desired value, make another slight adjustment in the relevant direction. **Be sure to change tips between volume setting adjustments and to pre-rinse each new tip.**

When the Pipetman F is readjusted to compensate for a particular solution, the original setting (ex. 0.2) can be noted on the self-adhesive labels supplied (Figure 4).

Figure 4

 *Liquids of high density, vapor pressure, or viscosity are best pipetted by the Gilson Microman® positive displacement pipette.*

5 - PRE-RINSING

When pipetting liquids having a viscosity and density different to water, for example organic solvents, a film of liquid is formed on the inside

wall of the pipette tip. This film can create an error. Since the film remains relatively constant in successive pipetting operations with the same tip, this error can be avoided by forming the film before transferring the first sample. This is done by aspirating a sample and dispensing it back into the same vessel. Since the film is already formed, all of the following samples will have better accuracy and repeatability.

This pre-rinsing operation should be repeated when a new tip is used.

6 - GILSON DIAMOND® TIPS

Gilson Diamond Tips are made to the highest specifications, strict quality control is maintained throughout the manufacturing process.


Diamond Tips are used to calibrate Pipetman F, therefore for optimum performance, you are strongly advised to use Gilson's Diamond Tips with your Pipetman F. Diamond Tips have the Gilson logo engraved on their collar, ensuring that you have a genuine Gilson product.

To ensure accuracy and precision, Gilson's Quality Assurance System focuses on the following critical parameters.

- Diamond Tips are made from pure polypropylene (virgin, metal-free, to avoid the possibility of contamination). They are available sterilized.
- Optimized shape (revised collar for optimum sealing, thin walls, and fine point), making them easier to mount, more flexible, with no vortexing, and improved precision.
- Diamond Tips are free from even microscopic defects, especially at the orifice. All surfaces are smooth and hydrophobic, thereby avoiding the excessive retention of liquids that causes poor accuracy and a lack of precision.

PIPETMAN® F

- Mold and cavity references are marked on the collar, ensuring the traceability for quality assurance purposes, batch numbers appear on all packages (bags and boxes).

- They form an air-tight seal with the tip-holder, preventing the leaks that cause poor accuracy and a lack of precision.
- They may be autoclaved at 121°C for 20 minutes at 0.1 MPa.

Diamond Filter Tips

Filter tips are used when sample-to-sample, sample-to-pipette, or sample-to-operator contamination must be avoided.

Because of their unique construction, Diamond Filter (DF) Tips may be autoclaved, without risk of damage to the filters. Gilson's sterilized Diamond Filter Tips are certified free of detectable RNase, DNase, DNA, RNA, and protease.

Diamond Tips to use for best results

Pipette	Diamond Tips	Type
F2 to F200	D200	Regular
F250 to F1000	D1000	Regular
Pipette	Diamond Filter Tips	Type
F2 to F25	DF30A DF30ST	Autoclavable Sterilized
F50 to F100	DF100A DF100ST	Autoclavable Sterilized
F200	DF200A DF200ST	Autoclavable Sterilized
F250 to F1000	DF1000A DF1000ST	Autoclavable Sterilized

7 - RECOMMENDATIONS

The recommendations below will ensure maximum performance from your Pipetman F.

- Make sure to operate the Pipetman F slowly and smoothly.
- The depth of immersion in the sample liquid should be the minimum necessary and should remain constant during aspiration.
- The Pipetman F should be held in the vertical position.
- Change the tip when a different liquid is to be aspirated, or when the Pipetman F is to be recalibrated.
- Change the tip if a droplet remains on the end of the tip from the previous pipetting operation.
- Each new tip should be pre-rinsed with the liquid to be pipetted.
- Liquid should never enter the Pipetman F tip-holder. To prevent this:
 - Press and release the push-button slowly and smoothly.
 - Never turn the pipette upside down.
 - Never lay the pipette on its side when there is liquid in the tip. Use a dedicated stand or hanger: Gilson Carrousel™, Trio™, or Single™.
- Never grease pistons or seals.

The pipette can be used between + 4 °C and + 40 °C but the specifications may vary according the temperature (see Chapter 2 for controlled conditions of use).

- When pipetting liquids with temperatures different from the ambient temperature, it is recommended to pre-rinse the tip several times before use.
- Do not pipette liquids with temperatures above 70 °C.
- After pipetting acids or corrosive solutions which emit vapours, it is recommended to disassemble the tip-holder and to rinse the piston and seal with distilled water.

8 - TROUBLESHOOTING

If there is a leak or the pipette is not accurate, check for the following:

- The connecting nut is loose.
Tighten the connecting nut.
- The tip-holder is cracked or scored.
Remove the tip ejector and inspect the tip-holder. For F2 to F20 models, if the tip-holder is damaged, the piston may also be damaged. Replace the damaged parts (see Chapter 9). When re-assembling the pipette, the connecting nut should be hand tightened.
- Chemical damage to the piston and seals.
Replace the piston and seals (see Chapter 9). Wash the inside of the tip-holder using distilled water.
- Improper reassembly.
Disassemble the pipette and reassemble the pipette correctly.

If there is liquid in the tip-holder, clean the pipette as follows:

- Remove the tip ejector. Unscrew the connecting nut and rinse the tip-holder, piston, seal and O-ring with distilled water. Dry these parts and reassemble the pipette.

 *If necessary, the tip-holder, connecting nut and the tip ejector can be autoclaved for 20 minutes at 121°C and at 0.1 MPa pressure. The seal and the O-ring should not be autoclaved.*

If there is an air bubble when the sample is aspirated:

- Eject the sample into its original vessel.
- Check that the tip is properly immersed in the sample liquid.
- Pipette the sample more slowly.

If the bubble appears a second time, replace the tip.

The outside of the Pipetman F can be cleaned using a cloth dampened with isopropanol.

If a problem continues after carrying out the above steps, send the pipette to your Gilson representative.

Before returning a pipette, please ensure that it is completely free of any chemical, biological, or radioactive contamination.

9 - REPLACEMENT PARTS


Description	F2 to F250	F30 to F100	F120 to F200	F250 to F1000
A Piston assembly	F123845	F144611	F123846	F123847
B Piston seal +	F144863	F144864	F144865	F144866
C O-ring				
D Connecting nut	F144723	F144723	F144723	F144723
E Ejector	F144766	F144767	F144768	F144769
F Tip-holder	F123353	F144602	F123305	F123371
G Calibration key	F123674	F123674	F123674	F123674

Contact your Gilson representative to order spare parts.

Replacing the piston will require re-calibration of the Pipetman F. We recommend that this work is carried out by your Gilson representative.

LT801118/J - ©2003 Gilson SAS All right reserved
English

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Gilson, Inc.