



Headoffice: C-4th&5th, 16, Deogyong-daero 1556beon-gil, Yeongtong-gu, Suwon-si, Gyeonggi-do, 16690, REPUBLIC OF KOREA
Manufacturing site: 74, Osongsaengmyeong 4-ro, Osong-eup, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, 28161, REPUBLIC OF KOREA

TEL : +82-31-300-0400 / FAX : +82-31-300-0499
<http://www.sdbiosensor.com>

Certificate of Analysis

| | | | |
|---------------------------|---|-----------------|----------------|
| Product Name | STANDARD™ F Influenza A/B FIA | | |
| Cat No. | 10INF20D | REF No. | FINF01G |
| Lot No. | FIN102009141 | MFG Date | 2020.09.14. |
| | | EXP Date | 2022.03.13. |
| Testing categories | Specification | | Result |
| Apperance Inspection | <ul style="list-style-type: none"> • Uniform Gap between upper and lower device • No scratch, stain and crack etc. • No contamination | | <i>Pass</i> |
| Sensitivity | In-house reference panels for the sensitivity test 1) All panels must be interpreted as positive. | | <i>Pass</i> |
| Specificity | In-house reference panels for the specificity test 1) All panels must be interpreted as negative. | | <i>Pass</i> |
| Precision | 1) Within-run precision : must be identified 100% agreement within-run test results 2) Between-run precision : must be identified 100% agreement between-run test results | | <i>Pass</i> |
| Packaging Inspection | Strip Devices (25strips) / Sterile swab (25ea) / Extraction buffer (25ea) / Filter cap(25ea)/Fixed volume dropper 300µl(25ea)/Positive control swab(1ea)/Negative control swab(1ea)/ insert (1ea) | | <i>Pass</i> |

According to SD Biosensor's quality assurance system, the above QC tests have been completed and passed the specification. We hereby approve this product release.

SD Biosensor, Inc

Geunkuk Song
 Quality Management Representative
 QA Division

Geunkuk Song
 Quality Management Representative
 QA Division