

STANDARD F Legionella Ag FIA




» Background

STANDARD F Legionella Ag FIA test system (Analyzer + Test device) detects *Legionella pneumophila* serogroup 1, 3, 5, 6 and 8 antigens via easily collectible urine sample. Without further sample processing, STANDARD F Legionella Ag FIA performs high sensitivity, and the test is not affected by Rheumatoid factor.

» Product Specification

STANDARD F Legionella Ag FIA	
Intended use	Qualitative detection of <i>L. pneumophila</i> serogroup 1, 3, 5, 6 and 8 antigens in urine specimens.
Specimen volume	100 µl
Specimen type	Urine
Measurement range	Quantitative result with Cut off index(COI)
Testing time	15 mins (Early detection available)
Storage condition	2 – 30°C / 36 - 86°F

» Test Procedure



» Test Performance

Reference	Sensitivity	Specificity
Fluorescent Immunoassay	94.9%(37/39)	100%(40/40)

» Ordering Information

Category	Product	Pack Size	CAT No.
Respiratory	STANDARD F Legionella Ag FIA	25 Tests	10LEG10D
Respiratory	Legionella Ag Control	Positive control x 10 Negative control x 10	10LEGC10